INTEGRA®

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Orotrakealt injektionsinstrument

## Product Information and Instructions

产品信息和说明 / 產品資訊及說明 / Informacije o proizvodu i upute / Informace o výrobku a pokyny / Produktinformation og -vejledninger / Productinformatie en instructies / Tuotetta koskevat tiedot ja ohjeet / Informations et instructions concernant le produit / Produktinformation und Gebrauchsanweisung / Informazioni sul prodotto e istruzioni per l'uso / 製品情報および使用説明 / Informacje o produkcie i instrukcje / Informações e instruções do produto / Informações e instruções acerca do produto / Информация об изделии и инструкции / Instrucciones e información sobre el producto / Produktinformation och instruktioner

Rx ONLY

# English

**Rx ONLY** Federal (USA) law restricts this device to sale by or on the order of a physician or practitioner

## Device Description / Intended Use

The Oro-Tracheal Injection Instrument is a syringe holder/injector intended to facilitate injections in the larynx or trachea from an oral approach. The Oro-Tracheal Injection Instrument is designed to be used with a special 27 gauge hypodermic needle (REF 1650050).

## Contraindications

Contraindications should not be used for anything other than their intended use.

## Warning

The Injection Instrument is provided non-sterile and must be cleaned and sterilized before the first use and any reuse.

DO NOT flash sterilize the oro-tracheal injector system. The instruments have not been validated for flash sterilization.

Discard instrument after suspected Creutzfeldt-Jakob Disease (CJD) exposure; the oro-tracheal injector system has not been validated to withstand the chemical and thermal exposures recommended to eradicate prions.

## Precautions

- USE ONLY SLIP-TIP SYRINGES. Syringes with locking type tips (i.e., Luer Lock) will not seat properly and the contents may leak.
- Do not use the syringe plunger to secure the syringe in place as the contents may be inadvertently discharged.
- All disposable needles and other sharps are to be immediately placed in an appropriate sharps disposal container at the point of use prior to transfer for decontamination

## Adverse Reactions

None known

## Pre-Use, Handling, And Inspection Of Instruments

- Read the Instructions for Use and keep them in a safe place.
- Use the product only in accordance with its intended use, see Indications.
- Use of an instrument for a task other than that for which it is intended could result in a damaged or broken instrument, or one which provides an unsatisfactory performance.
- Instruments should be handled and operated by personnel completely familiar with their use, assembly, and disassembly.
- Instruments must be rendered safe for handling, inspection, and assembly by wearing appropriate personal protection equipment (PPE) as promulgated by OSHA and AORN.
- Inappropriate use of instruments will lead to damage that is usually not repairable.
- Instruments must be thoroughly inspected upon receipt and prior to use to assure proper functioning. Failure to make a complete inspection to assure proper operation and function of instrument may result in unsatisfactory performance.
- The oro-tracheal injector system is supplied non-sterile and must be pre-cleaned, cleaned, visually examined, lubricated, and sterilized prior to use. Please see pre-cleaning, cleaning, visual examination, lubrication, and sterilization procedures below.
- Check screws on instruments after ultrasonic cleaning. Vibration from ultrasonic cleaning may cause them to loosen or fall out.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective.
- Replace any damaged components immediately with original spare parts.
- Store products in a dry, clean and safe place

## INSTRUCTIONS FOR USE

1. Select the needle cannula, curved or angled, best suited for the procedure to be performed.
2. Slide the threaded nut over the needle cannula and secure the cannula to the syringe handle/holder by turning the nut clockwise. Tighten the nut only enough to secure the cannula while allowing adjustment to the desired angle. Once the cannula is adjusted, tighten the nut until the cannula is secured to the handle/holder and will not rotate.
3. With the protective cap in place, screw the needle securely onto the threaded distal tip of the cannula by hand with a clockwise rotation. To prevent the possible leakage of fluids around the base of the needle hub, grasp the hub with hemostats and turn an additional turn to the right. Care should be exercised not to cross-thread the plastic-hubbed needle onto the metal threads of the cannula.
4. Insert a pre-filled 1.0 ml syringe with slip-tip into the handle holder. Make certain the slip-tip seats snugly by applying pressure on the syringe flange.

5. Once the syringe is securely seated and the cannula has been adjusted to the desired angle, remove the protective needle sheath. The injection instrument is ready for use. The syringe contents may be injected into the desired locations by exerting steady gentle pressure on the plunger in the same manner as with standard hypodermic syringes.
6. Following use, replace the protective needle sheath. Use hemostats to remove the needle by unscrewing it from the needle cannula in a counter-clockwise direction.
7. Dispose of the used syringe and needle safely in accordance with the health care facility's procedures.

Note: Needles are supplied sterile and are intended for single use only. Do not resterilize or reuse.

## Reprocessing Instructions for Reusable Instruments

When applying dry times to Integra LifeSciences instrument trays and their accessories, dry times outside the standard healthcare prevacuum parameters may be required.

### Drying Time:

Instruments Only and Aluminum/Stainless Steel Instrument Trays: 15 to 30 minutes and visual inspection for dryness.  
Polymer-based (plastic) Instrument Trays: Minimum drying time 40 minutes and visual inspection for dryness

<b>Reprocessing Instructions</b>	Cleaning and Sterilization instructions are subject to change.
<b>Warnings / Precautions</b>	To prevent stains, use distilled or demineralized water, and use a neutral detergent, to reprocess these instruments. Do not cold soak in glutaraldehyde, chlorine, or ammonium solutions, or dry sterilize, as damage to the instrument finish may occur. Hypodermic needles used with this injection instrument are intended for single use only, do not resterilize or reuse.
<b>Limitations</b>	After cleaning and sterilization, verify functionality prior to re-use.
<b>Point of Use</b>	This product is provided non-sterile and must be cleaned and sterilized before the first use and any reuse. Rinse instrument immediately after use in distilled or demineralized water. Keep instrument moist (for example, cover with a wet drape or saturate with tap water - do not use saline) after use and prior to decontamination, cleaning and sterilization to ensure adequate cleaning. Instruments with cleaning ports and/or suction tubes shall be flushed with neutral pH enzymatic detergent, followed by a tap water rinse.
<b>Containment and Transportation</b>	It is recommended that instruments are reprocessed as soon as is practical following use.
<b>Preparation for decontamination</b>	Fully disassemble modular designed instruments for effective cleaning. Remove any cap covering the cleaning port, if applicable. For instruments with lumens, and instruments equipped with cleaning ports, inject cleaning solution through the instrument with an irrigation syringe. When cleaning suction tubes with holes to control suction, place gloved finger over the control hole to flush through the tube. Open any articulated instruments before positioning in the soaking solution. Place instruments so they do not touch each other. <b>Note:</b> Flushing in running water is essential between decontamination and cleaning to prevent any risk of reaction between the two solutions.

Cleaning: Automated	<p>All instruments shall be thoroughly rinsed with tap water prior to placement in the automatic washer.</p> <p>Remove instruments and equipment from any sterilization trays before placing into washer baskets. Orient devices following recommendations of washer/disinfecting manufacturers.</p> <p>Use alkaline or neutral pH detergent recommended by washer/disinfecting or detergent manufacturers.</p> <p>These products have been validated for effective cleaning using an automatic washer/disinfecting cycle consisting of the following:</p>																			
	<table border="1"> <thead> <tr> <th>Phase</th> <th>Recirculation Time</th> <th>Water temperature</th> <th>Detergent type</th> </tr> </thead> <tbody> <tr> <td>Pre-wash 1</td> <td>3 minutes</td> <td>cold tap water</td> <td>N/A</td> </tr> <tr> <td>Enzyme Wash</td> <td>5 minutes</td> <td>hot tap water</td> <td>neutral pH enzymatic detergent</td> </tr> <tr> <td>Wash 1</td> <td>5 minutes</td> <td>66°C set point</td> <td>neutral pH detergent</td> </tr> <tr> <td>Rinse 1</td> <td>3 minutes</td> <td>hot tap water</td> <td>N/A</td> </tr> </tbody> </table> <p>Thoroughly examine instruments for any residual soil. A rigid MIS (minimally invasive surgery) rack should be used to clean rigid MIS instruments and instruments with lumens to maximize efficiency of cleaning. Suction ports and flush ports should be connected to the MIS Rack to allow a constant flow of solution inside hard-to-reach areas of the instruments.</p>	Phase	Recirculation Time	Water temperature	Detergent type	Pre-wash 1	3 minutes	cold tap water	N/A	Enzyme Wash	5 minutes	hot tap water	neutral pH enzymatic detergent	Wash 1	5 minutes	66°C set point	neutral pH detergent	Rinse 1	3 minutes	hot tap water
Phase	Recirculation Time	Water temperature	Detergent type																	
Pre-wash 1	3 minutes	cold tap water	N/A																	
Enzyme Wash	5 minutes	hot tap water	neutral pH enzymatic detergent																	
Wash 1	5 minutes	66°C set point	neutral pH detergent																	
Rinse 1	3 minutes	hot tap water	N/A																	
Cleaning: Manual	<p>Soak in lukewarm (less than 43°C), mild (pH 7.0 - 8.5), enzymatic detergent and deionized water for a minimum of 2 minutes. For instruments with lumens, and instruments equipped with cleaning ports, inject cleaning solution through the instrument with an irrigation syringe and allow to soak for a minimum of 2 minutes.</p> <p>Use a soft instrument brush to scrub instruments while submerged in cleaning solution to remove organic matter.</p> <p>Rinse with deionized water, then clean ultrasonically in a lukewarm (less than 43°C), mild (pH 7.0 - 8.5), detergent and deionized water for 10 minutes.</p> <p>Rinse thoroughly with deionized water, utilizing a syringe to thoroughly rinse cleaning solution from lumens and cleaning ports. Clean in this manner until no visible soil remains on the instrument.</p> <p>Dry with compressed air, or wipe dry with a lint-free cloth. Examine instruments for any staining or deterioration; remove from use as appropriate.</p> <p><b>Note:</b> When using an ultrasonic cleaner or a spray washing machine, follow the manufacturers recommendations, particularly with regard to articulated instruments and positioning of instruments.</p> <p>Following cleaning, lightly lubricate instruments with movable parts. Use a lubricant intended for sterilizable instruments such as a water-soluble instrument milk. Do not use silicone spray.</p>																			
Disinfection	<p><b>Note:</b> Do not cold soak in glutaraldehyde, chlorine, or ammonium solutions, or dry heat sterilize, as damage to the instrument finish may occur.</p>																			

Packaging	<p>A standard, sterilization wrap may be used. In the US, an FDA cleared surgical wrap is required.</p> <p>In sets: Instruments may be loaded into dedicated instrument trays or general purpose sterilization trays. Ensure that cutting edges are protected. Wrap trays using appropriate method.</p>																													
Sterilization	<p>Check the cleanliness and operation of the instrument. Clean again if debris is present and remove from use any damaged instrument. Close instruments with catches and racks on the first notch. Arrange the instruments in sterilization containers with perforations on the top and bottom, and on supports such as those used in microsurgery. Follow the appropriate cycle listed in the table below.</p> <p>All steam cycles have been validated in the wrapped configuration and can be sterilized wrapped or unwrapped. These devices have only been validated for steam sterilization methods.</p>																													
(Temperatures are minimum required, times are minimum required)																														
	<table border="1"> <thead> <tr> <th></th> <th colspan="2">Instruments Only</th> <th colspan="2">Instruments Only or Instruments in Instrument Trays</th> </tr> <tr> <th>Cycle</th> <th>Gravity</th> <th>Gravity</th> <th>Pre-vac</th> <th>Pre-vac (FR/WHO)</th> <th>Pre-vac (UK)</th> </tr> </thead> <tbody> <tr> <td>Temperature</td> <td>121°C</td> <td>132°C</td> <td>132°C</td> <td>134°C</td> <td>134°C</td> </tr> <tr> <td>Time</td> <td>30 minutes</td> <td>10 minutes</td> <td>4 minutes</td> <td>18 minutes</td> <td>3 minutes</td> </tr> <tr> <td>Drying</td> <td colspan="5"> <ul style="list-style-type: none"> <li>Instruments Only Aluminum/Stainless Steel Instruments Trays: 15 to 30 minutes and visual inspection for dryness.</li> <li>Poly-based (plastic) Instrument Trays: Minimum drying time 40 minutes and visual inspection for dryness.</li> </ul> </td> </tr> </tbody> </table>		Instruments Only		Instruments Only or Instruments in Instrument Trays		Cycle	Gravity	Gravity	Pre-vac	Pre-vac (FR/WHO)	Pre-vac (UK)	Temperature	121°C	132°C	132°C	134°C	134°C	Time	30 minutes	10 minutes	4 minutes	18 minutes	3 minutes	Drying	<ul style="list-style-type: none"> <li>Instruments Only Aluminum/Stainless Steel Instruments Trays: 15 to 30 minutes and visual inspection for dryness.</li> <li>Poly-based (plastic) Instrument Trays: Minimum drying time 40 minutes and visual inspection for dryness.</li> </ul>				
	Instruments Only		Instruments Only or Instruments in Instrument Trays																											
Cycle	Gravity	Gravity	Pre-vac	Pre-vac (FR/WHO)	Pre-vac (UK)																									
Temperature	121°C	132°C	132°C	134°C	134°C																									
Time	30 minutes	10 minutes	4 minutes	18 minutes	3 minutes																									
Drying	<ul style="list-style-type: none"> <li>Instruments Only Aluminum/Stainless Steel Instruments Trays: 15 to 30 minutes and visual inspection for dryness.</li> <li>Poly-based (plastic) Instrument Trays: Minimum drying time 40 minutes and visual inspection for dryness.</li> </ul>																													
Maintenance, Inspection and Testing	<p>Inspect components for any damage before and after each use. If damage is observed do not use the instrument until it is repaired. After cleaning and sterilization, verify functionality prior to re-use.</p>																													
Storage	<p>Store instruments in a clean, dry area.</p>																													
Additional Information	<p><b>Note:</b> Additional cleaning methods may be warranted, including presoaking in 3% hydrogen peroxide.</p>																													

Note: The instructions provided above have been validated by the manufacturer as being CAPABLE of preparing the product for re-use. It remains the responsibility of the processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the desired result. This normally requires validation and routine monitoring of the process.

#### CUSTOMER SERVICE INFORMATION

For further information regarding the use of this product or to report any problems, please contact Integra LifeSciences at 1-800-654-2873.

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Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries.
©2020 Integra LifeSciences Corporation. All Rights Reserved.
Integra LifeSciences Production Corporation