
Integra® Miltex®

Disposable Cervical Dilators

Directions for Use

Rx ONLY

REF 30-5001 Disposable Os Locator, Sterile

REF 30-5002 Disposable Canal Dilator, Sterile

REF 30-5003 Disposable Fundus Sound, Sterile

Description

The Os Locator

One molded band above the tactile handle and opposite the distal end of the device identifies the disposable Os Locator. The Os Locator may be used to find and gently dilate the external os.

The Canal Dilator

Two molded bands above the tactile handle identify the Canal Dilator from the Os Locator and the Fundus Sound. The Canal Dilator may be used to slowly dilate the length of the cervical canal.

The Fundus Sound

Three molded bands opposite the distal tip identify the Fundus Sound from the Os Locator and the Canal Dilator. The Fundus Sound is designed to pass through the entire length of the cervix, dilate the internal os and sound the fundus to determine uterine depth.

Indications

Cervical dilators are used in locating and dilating the external os, cervical canal and internal os. Cervical dilators are commonly used on post-menopausal patients with cervical stenosis.

Contraindications

Cervical dilators should not be used when the following conditions are present:

- Acute genital tract or pelvic infections.
- A pregnancy or the suspicion of a pregnancy.
- Any invasive cancer that is visible upon examination.
- Any non-compliant patient.

Precautions

- Check the packaging for damage which can compromise sterility.
- Visually inspect the dilators for any surface damage prior to use.
- Disposable cervical dilators are designed for single-patient use only.
- This procedure should only be performed by a trained medical professional.
- The use of excessive force should be avoided.
- Cervical dilators should not be used for prolonged periods of time.

Adverse Reactions

The following risks are possible during procedures requiring dilation:

- Laceration of the cervix
- Tears of the internal os
- Perforation of the fundus
- False passage at the internal os

Instructions

Dilation of the cervix may be required for many diagnostic and treatment procedures. This procedure should always be performed by a medical professional who is trained in the procedure and aware of all possible complications and risks that may be present.

1. A pelvic examination of the patient should be performed prior to dilation to determine conditions present and risks of dilation.
2. The insertion of a weighted speculum may assist with visualization of the cervix.
3. Uterine tenaculum forceps should be used to grasp the cervix. Stabilization of the tenaculum should be maintained during the dilation procedure to absorb resistance and avoid trauma to surrounding anatomical structures.
4. The Os Locator is used to gently dilate the external os to allow access to the cervical canal. Holding the dilator like a pen, gently insert the dilator into the external os. Slight, slow movements can be used to break up connective tissues that may exist. Proceeding slowly allows the os to gently stretch to accommodate the dilator.
5. The Canal dilator is often needed next to continue the dilation of the cervical canal to gain access to the internal os. The canal may become atrophic and rigid in post menopausal patients. Proceeding slowly using the technique referenced above, allows the canal to gently stretch over the dilator.
6. The Fundus Sound is utilized to complete the dilation of the internal os and palpate the depth of the uterus. Continuing to hold the dilator like a pen, dilation of the internal os should also proceed slowly to allow the tissues to gently stretch over the dilator. Care should be exercised to minimize the risk of perforation of the uterus or tearing of the internal os.
7. If additional dilation is required, the use of progressively larger diameter dilators can be utilized following the same process allowing for slow dilation of the cervix to accommodate instrumentation needed for the procedure.
8. If resistance is met to any size dilator, it is advisable to leave it in the internal os for several seconds to allow the tissues to expand around it.

Symbols

Symbols may be used on some package labeling for easy identification.

Manufacturer¹

Catalog number

Lot Number

See instructions for use - Caution

Consult IFU and website address

Single use only

Expiration date

Method of sterilization: ethylene oxide

Federal (USA) law restricts this device to sale by or on the order of a physician or practitioner

¹ Company responsible for a device marketed under its own name regardless of whether "manufactured for" or "manufactured by" the company.

Manufactured for
Integra LifeSciences Production Corporation
11 Cabot Boulevard
Mansfield, MA 02048 USA
integralife.com

Manufacturer
PANPAC Medical Corporation
2~6F, No.58, Nanxing Rd,
Shi-chih Dist., New Taipei City,
221 Taiwan, R.O.C

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