Integra®

OSV II® Lumbar Valve Systems

Systèmes de Valve Lombaire OSV II®

OSV II® Lumbale Ventilsysteme

Sistema Lombare a Valvola OSV II®

Sistemas de Válvula Lombar OSV II®

Sistemas de Válvula Lumbar OSV II®

OSV II® Lumbale Klepsystemen

Instructions for Use

Mode d'Emploi

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MANUFACTURER:


Integra LifeSciences Switzerland Sarl

Rue Girardet 29 (2nd Floor)

Le Locle CH-2400, Switzerland

Figure 1 - Schematic Drawing – Integra® OSV II® Lumbar Valve Mechanism

- Figure 1 - Dessin schématique - Mécanisme de la valve lombaire Integra® OSV II®
- Abb. 1 - Schematische Darstellung - Integra® OSV II® Lumbar-ventilmechanismus
- Figura 1 - Disegno schematico del meccanismo della valvola lombare Integra® OSV II®
- Figura 1 - Desenho esquemático - Mecanismo da válvula Lombar Integra® OSV II®
- Figura 1 - Esquema del mecanismo de válvula Lumbar Integra® OSV II®
- Afb. 1 - Schematische tekening – Mechanisme van de Integra® OSV II® Lumbale klepsystemen


1. Inlet/Entrée/Einlaß/Entrata/Entrada/Entrada/Ingang
2. Outlet/Sortie/Auslaß/Uscita/Saída/Salida/Uitgang
3. Case/Boîtier/Gehäuse/Involucro/Involucro/Carcasa/Behuizing
4. Diaphragm/Diaphragme/Membran/Diaframma/Diafragma/Diafragma/Diafragma
5. Seat/Siège/Rubining/Alloggiamento/Assento/Asiento/Zitting
6. Notched pin/Pointeau/Rubinstift mit Aussparung/Perno/Pino dentado/Pasador dentado/Gekerfde pin

Figure 2 - Schematic Drawing - Example of Integra® OSV II® Lumbar Valve Flow/Pressure Curve


- Figure 2 - Schéma: Exemple de courbe pression/débit de la valve lombaire Integra® OSV II®
- Abb. 2 - Schematische Darstellung Integra® OSV II® Lumbar ventil fluss/druk-kurve
- Figura 2 - Disegno schematico della curva flusso/pressione della valvola lombare Integra® OSV II®
- Figura 2 - Desenho esquemático - Exemplo da curva de pressão/fluxo da válvula Lombar Integra® OSV II®
- Figura 2 - Ejemplo de curva de variación de la tasa de flujo/gradiente de presión para la válvula Lumbar Integra® OSV II®
- Afb. 2 - Schematische tekening- Voorbeeld van de stroom/druk curve van de Integra® OSV II® Lumbar klepsystemen

Differential pressure / Pression différentielle / Differentialdruck / Presione differenziale
 Pressão diferencial / Presión diferencial / Verschilddruk (mmH₂O)



- ① Stage 1/Phase 1/Stufe 1/Stadio 1/Fase 1/Etapa 1/Fase 1 : **Low differential pressure**/Basse pression différentielle/
 Niedriger Differenzdruck/Bassa pressione differenziale/Pressão diferencial baixa/Gradiente de presión bajo/Lage verschilddruk
- ② Stage 2/Phase 2/Stufe 2/Stadio 2/Fase 2/Etapa 2/Fase 2 : **Flow regulation**/Régulation de débit/Flußregulierung/
 Regolazione del flusso/Regulação do fluxo/Regulación de la tasa de flujo/Stroomregulering
- ③ Stage 3/Phase 3/Stufe 3/Stadio 3/Fase 3/Etapa 3/Fase 3 : **Pressure relief mode**/Limitation de pression/
 Sicherheitsstufe/Limitazione della prezione/Alívio de pressão/Fase de seguridad/Drukbeperking

- Figure 3** - Patency Testing Procedure
Figure 3 - Procédure de vérification de la perméabilité
Abb. 3 - Überprüfung der Ventildurchlässigkeit
Figura 3 - Procedura per il test di pervietà
Figura 3 - Procedimento do teste de desobstrução
Figura 3 - Procedimiento de comprobación de ausencia de obstrucción
Afb. 3 - Testprocedure voor de dooorlaatbaarheid

- Figure 4** - X-Ray image
Figure 4 - Image radiographique
Abb. 4 - Röntgenbild
Figura 4 - Immagine a raggi X
Figura 4 - Reprodução radiológica
Figura 4 - Imagen de rayos X
Afb. 4 - Röntgenfoto

DESCRIPTION

The Integra® OSV II® Lumbar hydrocephalus Valve System is an implantable device for controlled cerebrospinal fluid (CSF) drainage from the lumbar subarachnoid region to the peritoneal cavity. Unlike conventional valves, the OSV II Lumbar Valve System is a variable resistance valve which maintains drainage rate close to the CSF secretion rate (20 ml/hr) within the physiological range of intracranial pressure (ICP). The mechanism incorporates a safety pressure relief mode to prevent accidental intracranial hypertension.

Caution: The Integra OSV II Lumbar Valve may appear similar to the Integra® Flow Regulating Valve Low Flow, Lumbar, however, the Integra Flow Regulating Valve Low Flow, Lumbar features a black, radiopaque flow direction arrow to differentiate both products: visually during implantation and by X-ray after implantation (Figure 4).

WARNING:

The Integra OSV II Lumbar Valve should be used in patients who need a drainage rate close to the CSF secretion rate (18-30 ml/hr). For patients who need a lower drainage rate (8-17 ml/hr), the Integra Flow Regulating Valve Low Flow, Lumbar should be used.

Principles of Operation

The variable flow restrictor (Figure 1) consists of a silicone elastomer diaphragm (4), synthetic ruby seat (5), and notched pin (6). The diaphragm reacts to differential pressure (DP) variations. The seat is inserted into the center of the diaphragm. As differential pressure varies, the clearance between the seat and the pin increases or decreases depending on seat movement along the pin. Changes in the seat and pin clearance vary flow rates. A typical flow/pressure curve is shown in Figure 2.

Note - Differential pressure/flow characteristics for each valve are verified to be within specification at time of manufacture.

The three stages of operation of the OSV II Lumbar Valve System are defined below:

Stage I Low Differential Pressure

This stage begins when the flow rate through the valve reaches 5ml/hr (DP will be between 30 and 80 mm H₂O). The valve remains in Stage I with CSF flow rates up to 18 ml/hr (DP will be between 40 and 120 mm H₂O).

Stage II Flow Regulation

When the DP increases, the valve operates as a variable resistance flow regulator. At DP ranges between 120 and 300 mm H₂O, the valve restricts flow between 18 and 30 ml/hr.

Stage III Pressure Relief Mode

Should the intraventricular pressure (IVP) elevate abruptly, the shunt operates in a rapid flow rate mode to facilitate IVP normalization. The valve then reverts to Stage II or I, depending upon conditions.

Configurations

The system is available in two configurations: valve with or without antechamber (see product label). The system consists of a valve with pre-attached, proximal tubing and pre-attached distal catheter, one F5 lumbar catheter, one guidewire, one 14G Tuohy needle, one suture clamp, one step-down connector, one luer connector and one straight connector.

Materials

Materials in contact with body tissues or body fluids are silicone elastomer with or without barium sulfate, polysulfone, synthetic ruby, silicone adhesive, polypropylene, epoxy, stainless steel, polycarbonate and PTFE.

INDICATIONS

The OSV II Lumbar Valve System is an implantable system used in the treatment of patients with communicating hydrocephalus to shunt CSF from the lumbar subarachnoid region to the peritoneal cavity.

CONTRAINDICATIONS

The OSV II Lumbar Valve System should not be used in patients with non-communicating hydrocephalus.

This valve system should not be implanted when an infection along the shunt pathway (e.g. meningitis, ventriculitis, peritonitis) is suspected. It is advisable to postpone valve implantation if an infection is present anywhere in the body (septicemia, bacteremia).

The OSV II Lumbar Valve System should not be implanted in patients with untreated choroid plexus tumors. Such tumors produce CSF at rates in excess of the specification of the flow regulation Stage II, and the OSV II Lumbar Valve System would underdrain under these conditions.

SIDE EFFECTS

- In addition to the risks associated with shunt system or component implantation, major complications include mechanical failure, such as tubing fracture, material deterioration with time, or shunt pathway obstruction, infection, immune system reaction to materials and CSF leakage along the shunt pathway. System disconnection can lead to catheter migration into the subarachnoid space or peritoneal cavity.
- Intracranial hypertension syndrome is a specific result of obstruction.
- Overdrainage may result in subdural hematomas, craniostenosis, intracranial hypotension syndrome, sunken fontanelles (in infants) or acquired tonsillar herniation.
- Shunt system failure should be suspected if any of these is noted: Continuing symptoms of hydrocephalus; CSF leakage; redness, tenderness or erosion of the skin along the shunt pathway. These complications require prompt replacement of the shunt system or affected component(s).
- Low grade bacterial colonization may cause recurrent fevers, anemias, splenomegaly, shunt nephritis or pulmonary hypertension.

- Possible complications of lumboperitoneal shunting include arachnoiditis, scoliosis, hyperlordosis, and back pain.
- Perforation of an abdominal viscus is a rare complication.

PRECAUTIONS

- Closely observe patients with an implanted valve system for shunt failure symptoms (see "Side Effects").
- During implantation, have duplicates of each component available.
- Handle valve tubing carefully. Cover forceps and clamp jaws with silicone elastomer tubing.
- Do not lubricate valves. They will be lubricated adequately by CSF.
- Silicone elastomer requires special handling to prevent contamination. Avoid contact with glove talc or skin oils.
- Other than trimming catheters, modifications are not recommended.
- The kits contain the non-implanted accessories used to aid in implantation. **Do not** bring these accessories into the MR environment as serious injury can occur.
- **Note** - 10 cm of the Lumbar catheter will add, at 10 ml/hr, an extra resistance of 3.1 mmH₂O (hydrodynamic resistance of the Lumbar catheter: 0.031 mm H₂O per ml/hr and per cm).

Caution: Do not overpressurize the valve system. Excessive flushing pressures may lead to valve damage.

Caution: Do not use in conjunction with other devices, such as antisiphon devices. Addition of such components change the flow/pressure characteristics of the valve.

Caution: This product is for single use only. Re-sterilization and re-use may affect its performances and characteristics and may present a contamination risk (such as Creutzfeldt-Jakob disease).

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or practitioner.

Storage and Transport: The valve systems in Integra packaging may be stored at room temperature. They will withstand the conditions of normal transport. *This product is sterilized with ethylene oxide. Do not use if the package is open or damaged. Use the device prior to the «Use by» date on the package label.*

Caution - Do not resterilize. Integra will not be liable for any or all damages including, but not limited to, direct, indirect, incidental, consequential or punitive damages resulting from or related to resterilization.

OSV II LUMBAR VALVE PRESSURE / FLOW TEST

Each OSV II Lumbar Valve System is tested at the end of the manufacturing process. Integra does not recommend testing the valve in the operating room before implantation due to risk of possible contamination and because the conditions required for this testing are difficult to meet in a clinical environment.

PREIMPLANT TESTING FOR PATENCY

Should verification of valve patency be desired prior to implantation, place the inlet of the valve system in a container of sterile apyrogenic saline solution. Connect a syringe, with the Luer connector attached to the drainage catheter, and gently draw the solution through the system to displace air (**Figure 3**). Remove syringe and Luer connector. Holding the valve's inlet end in the sterile solution, allow the drainage catheter to hang freely (**Figure 3**). The valve is patent if solution drips from the drainage catheter. Time between each drop will be approximately 4-6 seconds.

Note: This procedure only tests the valve patency. It does not allow verification of the 3 different stages of operation of the OSV II Lumbar Valve System.

RECOMMENDED IMPLANTATION PROCEDURE

Before implantation, the valve serial number should be recorded in the patient's file.

Note: After implantation, complete the patient identification card for the patient or family.

Clinical Judgment

Incision size, size and type of nonabsorbable suture material, etc. should be selected according to the surgeon's experience and preference. Procedure may be varied in accordance with the surgeon's clinical judgment.

Modifications

The lumbar and drainage catheters can be trimmed to the desired length during the implantation procedure. The lumbar catheter should be trimmed only at its distal end since the proximal end contains drainage holes necessary for proper function. **Other modifications are not recommended.**

Preparation

1. Administer appropriate anesthesia.
2. Place the patient in the lateral position, left side up. Instruct the patient to flex the right (lower) leg and to extend the left (upper) leg.
3. Prepare and drape the lower back, hip and abdomen.
4. Make an incision above the iliac crest.

Placement of the lumbar catheter

Caution: Do not use bacitracin solution to flush the catheter and guidewire as it may impair the insertion and removal of the guidewire. Use only sterile saline solution to flush the catheter and guidewire.

1. At the physician's option, make a 2 cm incision prior to Tuohy needle insertion.
2. Connect a syringe to the Luer-lock hub of the guidewire dispenser and flush the guidewire with sterile saline solution.
3. Insert the Luer connector into the catheter. Connect a syringe to the adapter and flush the catheter with sterile saline solution.
4. Insert the guidewire, PTFE-coated tip first, into the catheter.
5. Insert the Tuohy needle, with obturator in place, into the subarachnoid space. Insert the needle bevel up for cephalad placement of the catheter; bevel down for caudad placement.
6. To verify that the lumen of the needle has entered subarachnoid space, withdraw the obturator slowly until CSF appears. If it does not, reinsert the obturator and reposition the needle.
7. Verify that the needle is correctly positioned. Then, remove the obturator and insert the catheter and guidewire through the needle to the desired lumbar level.
8. Remove the Luer connector and withdraw the needle over the end of the catheter.
Caution: To avoid damage to the catheter, do not withdraw the catheter after it has been inserted into the needle. If the catheter must be removed, withdraw the needle and catheter simultaneously.
9. Remove the guidewire. **Caution:** To avoid damage to the catheter, do not withdraw the guidewire without first removing the Tuohy needle.
10. **Connect a syringe with the provided Luer connector to the lumbar catheter and aspirate 2 to 3 ml of CSF to eliminate possible debris from the CSF.**
11. Clamp the lumbar catheter at incision level.
12. Anchor the lumbar catheter to fascia using the suture clamp.
13. Using a tunneler or other suitable instrument, pass the lumbar catheter subcutaneously from the back to the previous incision. **Do not reduce catheter length at this time.**

Placement of the valve

1. Fill the valve system with a sterile apyrogenic saline solution by gentle aspiration through the distal catheter.
2. Using blunt dissection, create a small pocket under the skin to avoid unnecessary pull strength on the drainage catheter when positioning the shunt.
3. Trim the free (output) end of the lumbar catheter to the desired length. Connect the lumbar catheter to the tubing of the valve using the stepdown connector, place the small end in lumbar catheter. For easy connection, insert the guidewire through the stepdown connector and a short distance into the distal end of the catheter. Slide the catheter onto the connector and remove the guidewire.
4. Secure all connections with non absorbable sutures, and tie catheter to underlying fascia.
5. To check valve patency, release the clamp on the lumbar catheter and remove air from the system by allowing CSF to flow through the catheter.
6. **Caution:** Prior to the introduction of the drainage catheter into the abdominal cavity, it is advisable to confirm the patency of the system by observing the formation of CSF droplets at the end of the drainage catheter.
Notes:
 - The CSF droplets formation rate may be lower than that observed with conventional DP shunts because of the flow restriction mode of the OSV II Lumbar Valve System.
 - In children, insert sufficient drainage catheter length into the abdominal cavity to allow for growth.
7. Introduce the drainage catheter in the peritoneal cavity using a trocar or standard laparotomy technique consistent with the surgeon's experience.
8. Close the superficial fascia and skin in a routine manner.

Note - X-ray the complete system just after implantation for future reference to determine whether system components have shifted.

Caution: When sampling or infusing solution using the antechamber, a 25G or smaller beveled needles should be used. Care must be taken during insertion and removal of the needle to avoid tearing the silicone material. It is also advised to avoid puncturing several times at the same place. Taking into account the above recommendations, the antechamber is qualified to sustain up to 40 punctures without leaking.

Caution: Do not overpressurize the valve system. Excessive flushing pressures may lead to valve damage.

MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that the implanted OSV II Valve systems components are MR CONDITIONAL. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3.0 T or less
- Maximum spatial field gradient of 3500 gauss/cm (35 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Mode)

Do not take non-implanted accessories used to aid in introduction of the shunt device such as stylets and catheter guide into the MR environment. They have not been evaluated for safety in the MR environment.

Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than 0.4 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the OSV II Valve system extends approximately 10 mm from the device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

The flow characteristics of the shunt are not expected to be affected by MR exposure because the valve mechanism does not have any metallic components.

HOW SUPPLIED

Integra OSV II Lumbar Valve Systems are supplied sterile and non-pyrogenic in a double-wrap packaging.

Disposal

After patient use, the system must be handled as biohazardous material and disposed of in accordance with applicable Federal, State, local, or international environmental requirements following facility protocols.

PRODUCT INFORMATION DISCLOSURE

INTEGRA HAS EXERCISED REASONABLE CARE IN THE CHOICE OF MATERIALS AND MANUFACTURE OF THIS PRODUCT. INTEGRA EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. INTEGRA SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THIS PRODUCT. INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS. INTEGRA INTENDS THAT THIS DEVICE SHOULD BE USED ONLY BY PHYSICIANS WITH EDUCATIONAL AND TRAINING BACKGROUND ENABLING THE PROPER USE OF THE DEVICE.

Returned Goods Policy

Authorization from customer service must be obtained prior to returning the product.

Sterile products must be returned in unopened packages with manufacturer's seals intact to be accepted for replacement or credit unless returned due to a complaint or product defect. Determination of a product defect will be made by Integra. Credit will be issued for goods returned prior to 90 days from ship date and may be subject to a restocking charge of up to 20%. This assumes that the product returned is not damaged and can be verified to have not been used or opened.

Product Order Information

All products can be ordered through your Integra NeuroSpecialist or customer service representative or by contacting:

Integra LifeSciences Corporation
1100 Campus Road
Princeton, NJ 08540, U.S.A.
Telephone: 1-800-654-2873
Outside the US: +1-609-275-0500
Fax: +1-609-275-5363

Service and Repair
For service, repair or replacements, contact your
local Integra representative.

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Integra, the Integra logo, and OSV II are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries.
© 2022 Integra LifeSciences Corporation. All Rights Reserved.