Integra®

OSV II® Valve Unit

OSV II® Valve System

OSV II® Low Pro Valve Unit

OSV II® Low Pro Valve System


Instructions for Use
Mode d'Emploi
Gebrauchsanleitung
Istruzioni per l'Uso
Instruções de utilização
Instrucciones de Uso
Gebruiksaanwijzing

MANUFACTURER:
Integra LifeSciences Switzerland Sarl
Rue Girardet 29 (2nd Floor)
Le Locle CH-2400, Switzerland

Figure 1 - Schematic Drawing - Integra® OSV II® Valve Mechanism

- Figure 1 - Dessin schématique - Mécanisme de la valve Integra® OSV II®
- Abb. 1 - Schematische Darstellung - Integra® OSV II® Ventilmechanismus
- Figura 1 - Disegno schematico del meccanismo della valvola Integra® OSV II®
- Figura 1 - Desenho esquemático - Mecanismo da válvula Integra® OSV II®
- Figura 1 - Esquema del mecanismo de válvula Integra® OSV II®
- Afb. 1 - Schematische tekening - Mechanisme van de Integra® OSV II® klep


1. Inlet/Entrée/Einlaß/Entrata/Entrada/Ingang
2. Outlet/Sortie/Auslaß/Uscita/Saída/Salida/Uitgang
3. Case/Boîtier/Gehäuse/Involucro/Invlucro/Carcasa/Behuizing
4. Diaphragm/Diaphragme/Membran/Diaframma/Diafragma/Diafragma/Diafragma
5. Seat/Siège/Rubining/Alloggiamento/Assento/Asiento/Zitting
6. Notched pin/Pointeau/Rubinstift mit Aussparung/Perno/Pino dentado/Pasador dentado/Gekerfde pin

Figure 2 - Schematic Drawing - Example of Integra® OSV II® valve Flow/Pressure Curve


- Figure 2 - Schéma: Exemple de courbe pression/débit de la valve Integra® OSV II®
- Abb. 2 - Schematische Darstellung Integra® OSV II® ventil fluss/druck-kurve
- Figura 2 - Disegno schematico della curva flusso/pressione della valvola Integra® OSV II®
- Figura 2 - Desenho esquemático - Exemplo da curva de pressão/fluxo da válvula Integra® OSV II®
- Figura 2 - Ejemplo de curva de variación de la tasa de flujo/gradiente de presión para la válvula Integra® OSV II®
- Afb. 2 - Schematische tekening- Voorbeeld van de stroom/druk curve van de Integra® OSV II® klep

Differential pressure / Pression différentielle / Differentialdruck / Presione differenziale / Pressão diferencial / Presión diferencial / Verschuldruk (mmH₂O)


1. **Stage 1/Phase 1/Stufe 1/Stadio 1/Fase 1/Etapa 1/Fase 1 : Low differential pressure/Basse pression différentielle/ Niedriger Differenzdruck/Bassa pressione differenziale/Pressão diferencial baixa/Gradiente de presión bajo/Lage verschuldruk**
2. **Stage 2/Phase 2/Stufe 2/Stadio 2/Fase 2/Etapa 2/Fase 2 : Flow regulation/Régulation de débit/Flußregulierung/Regolazione del flusso/Regulação do fluxo/Regulación de la tasa de flujo/Stroomregulering**
3. **Stage 3/Phase 3/Stufe 3/Stadio 3/Fase 3/Etapa 3/Fase 3 : Pressure relief mode/Limitation de pression/ Sicherheitsstufe/Limitazione della prezione/Alivio de pressão/Fase de seguridad/Drukbeperking**


- Figure 3** - Patency Testing Procedure
 Figure 3 - Procédure de vérification de la perméabilité
 Abb. 3 - Überprüfung der Ventildurchlässigkeit
 Figura 3 - Procedura per il test di pervietà
 Figura 3 - Procedimento do teste de desobstrução
 Figura 3 - Procedimiento de comprobación de ausencia de obstrucción
 Afb. 3 - Testprocedure voor de dooorlaatbaarheid

- Figure 4** - Ventricular Catheter Introducer
 Figure 4 - Introducteur de cathéter ventriculaire
 Abb. 4 - Zweiteilige Einführhilfen für integrierte Ventrikelkatheter
 Figura 4 - Introduttore del catetere ventricolare
 Figura 4 - Introdudor do cateter ventricular
 Figura 4 - Introducitor del Catéter Ventricular
 Afb. 4 - Ventriculaire Katheteregeleider


- Figure 5** - Antechamber occlusion for pumping
 Figure 5 - Utilisation de l'antichambre pour pompage du système
 Abb. 5 - Vorkammer okklusion für Pumpen
 Figura 5 - Utilizzazione dell'anticamera per pompaggio di sistema
 Figura 5 - Utilização da antecâmara para bombeamento do sistema
 Figura 5 - Utilisación della antecámara por bombeo del sistema
 Afb. 5 - Occlusie voorruimte om te pompen


- ① Proximal side occlusion/Obstruction côté proximal/Okklusion seide proximal/Occlusione estremità prossimale/
 Oclusão do lado proximal/Occlusion lado proximal/Occlusie proximale kant
- ② Distal side occlusion/Obstruction côté distal/Okklusion seide distal/Occlusione estremità distale/
 Oclusão do lado distal/Occlusion lado distal/Occlusie distal kant

- Figure 6 - Schematic Drawing - Top view
 Figure 6 - Dessin schématique - Vue de dessus
 Abb. 6 - Schematische Darstellung - Ansicht von oben
 Figura 6 - Disegno schematico - Vista dall'alto
 Figura 6 - Desenho esquemático - Vista de cima
 Figura 6 - Esquema - Vista superior
 Afbeelding 6 - Schematische tekening - Bovenanzicht

Integra® OSV II® Low Pro Valve


Integra® OSV II® Valve


- Figure 7 - Schematic Drawings - Side view
 Figure 7 - Dessin schématique - Vue de profil
 Abb. 7 - Schematische Darstellung - Profilsicht
 Figura 7 - Disegno schematico - Vista laterale
 Figura 7 - Desenho esquemático - Vista lateral
 Figura 7 - Esquema - Vista lateral
 Afbeelding 7 - Schematische tekening - zijaanzicht

Integra® OSV II® Low Pro Valve



Integra® OSV II® Valve

- ① Antechamber/ Vorkammer/ anticamera/ antecâmara/ antecámara/ voorruimte
 ② Valve mechanism/ Mécanisme de la valve/ Ventilmechanismus / Meccanismo della valvola/ Mecanismo da válvula/
 Mecanismo de la válvula/ klepmechanisme

- Figure 8 - X-Ray image
 Figure 8 - Image radiographique
 Abb. 8 - Röntgenbild
 Figura 8 - Immagine a raggi X / Reprodução radiológica / Imagen de rayos X
 Afbeelding 8 - Röntgenfoto

DESCRIPTION

The OSV II® Hydrocephalus Valve System (OSVII, OSVII Low Pro) is an implantable device for controlled cerebrospinal fluid (CSF) drainage from the ventricles to the peritoneal cavity or other appropriate drainage site such as the heart's right atrium. Unlike conventional valves, the OSV II Valve System is a variable resistance valve which maintains a drainage rate close to the CSF secretion rate (20 ml/hr) within the physiological range of intracranial pressure (ICP). The mechanism incorporates a safety pressure relief mode to prevent accidental intracranial hypertension.

Caution: The Integra® OSV II® Valve may appear similar to the Integra® Flow Regulating Valve Low Flow, however, the Integra Flow Regulating Valve Low Flow features a black radiopaque flow direction arrow to differentiate both products, visually during implantation, and by X-ray, after implantation (See OSV II Valve X-ray image in Figure 8).

WARNING:

The Integra OSV II Valve System should be used in patients who need a drainage rate close to the CSF secretion rate (18-30 ml/hr). For patients who need a lower drainage rate (8-17 ml/hr), the Integra Flow Regulating Valve Low Flow should be used.

Principles of Operation

The variable flow restrictor (Figure 1) consists of a silicone elastomer diaphragm (4), synthetic ruby seat (5), and notched pin (6). The diaphragm reacts to differential pressure (DP) variations. The seat is inserted into the center of the diaphragm. As differential pressure varies, the clearance between the seat and the pin increases or decreases depending on seat movement along the pin. Changes in the seat and pin clearance vary flow rates. A typical flow/pressure curve is shown in Figure 2.

Note - Differential pressure/flow characteristics for each valve are verified to be within specification at time of manufacture.

The three stages of operation of the OSV II Valve System are defined below:

Stage I Low Differential Pressure

This stage begins when the flow rate through the valve reaches 5ml/hr (DP will be between 30 and 80 mm H₂O). The valve remains in Stage I with CSF flow rates up to 18 ml/hr (DP will be between 40 and 120 mm H₂O).

Stage II Flow Regulation

When the DP increases, the valve operates as a variable resistance flow regulator. At DP ranges between 120 and 300 mm H₂O, the valve restricts flow between 18 and 30 ml/hr.

Stage III Pressure Relief Mode

Should the intraventricular pressure (IVP) elevate abruptly, the shunt operates in a rapid flow rate mode to facilitate IVP normalization. The valve then reverts to Stage II or I, depending upon conditions.

Configurations

The system is available in different configurations (see product label). A system generally consists of a ventricular (proximal) catheter, a valve unit, and a drainage (distal) catheter. Depending on the configuration, the ventricular catheter is fully radiopaque with radiopaque length markings at 2cm intervals from the tip or with a barium stripe and non radiopaque length markers every 1cm. Some configurations contain an antechamber (standard or low profile "Low Pro", see Figures 6, 7, 8) or a burr hole cap. Connectors are designed for use with 1.1 to 1.4 mm inner diameter silicone elastomer tubings.

Implanted Materials

The implanted product may contain one or more of the following: Silicone elastomer with or without barium sulfate, silicone ink, polysulfone, synthetic ruby, silicone adhesive, tantalum, polypropylene and epoxy.

INDICATIONS

The OSV II Valve System is an implantable system used in the treatment of patients with hydrocephalus, to shunt CSF from the ventricles to the peritoneal cavity or other appropriate drainage site such as the heart's right atrium.

CONTRAINDICATIONS

This valve system should not be implanted when an infection along the shunt pathway (e.g. meningitis, ventriculitis, peritonitis) is suspected. It is advisable to postpone valve implantation if an infection is present anywhere in the body (septicemia, bacteremia). Atrial shunting is not advised for patients with congenital heart disease or other serious cardiopulmonary abnormalities.

The OSV II Valve System should not be implanted in patients with untreated choroid plexus tumors. Such tumors produce CSF at rates in excess of the specification of the flow regulation Stage II, and the OSV II Valve System would underdrain under these conditions.

The OSV II Valve System should not be used for drainage of extraventricular fluid collections such as hygromas or cysts; such conditions are typically treated with very low differential pressure valves.

SIDE EFFECTS

- In addition to the risks associated with shunt system or component implantation, major complications include mechanical failure, such as tubing fracture, material deterioration with time, or shunt pathway obstruction, infection, immune system reaction to materials and CSF leakage along the shunt pathway. System disconnection can lead to catheter migration into the atrium, peritoneum or lateral ventricles.
- Ventricular shunting has been associated with epilepsy, and shunt dependency. Intracranial hypertension syndrome is a specific result of obstruction. Overdrainage may result in subdural hematomas, craniostenosis, intracranial hypotension syndrome or sunken fontanelles (in infants).
- Shunt system failure should be suspected if any of these is noted: Continuing symptoms of hydrocephalus; CSF leakage; redness, tenderness or erosion of the skin along the shunt pathway. These complications require prompt replacement of the shunt system or affected component(s).
- Low grade bacterial colonization may cause recurrent fevers, anemias, splenomegaly, shunt nephritis or pulmonary hypertension.
- With atrial catheterization, the possibility of clot formation around the atrial catheter tip could lead to pulmonary tree embolization.
- Perforation of an abdominal viscus is a rare complication of ventriculoperitoneal shunts.

PRECAUTIONS

- Closely observe patients with an implanted valve system for shunt failure symptoms (see "Side Effects").
- During implantation, have duplicates of each component available.
- Handle valve tubing carefully. Cover forceps and clamp jaws with silicone elastomer tubing.
- Do not lubricate valves. They will be lubricated adequately by CSF.
- Silicone elastomer requires special handling to prevent contamination. Avoid contact with glove talc or skin oils.
- The drainage catheter may be trimmed to the desired length during implantation. Other modifications are not recommended.
- To minimize the possibility of underdrainage in patients confined to a supine position (i.e., patients under 3 months of age and bedridden patients), these patients should be elevated at a 45° angle for several hours a day (using a pillow, etc.) until they can sit or stand independently.
- Some kits contain the non-implanted accessories used to aide in implantation. **Do not** bring these accessories into the MR environment as serious injury may occur.

Caution: Do not overpressurize the valve system. Excessive flushing pressures may lead to valve damage.

Caution: Do not use in conjunction with other devices, such as antisiphon devices or distal slit valves. Addition of such components change the flow/pressure characteristics of the valve.

Caution: This product is for single use only. Re-sterilization and re-use may affect its performances and characteristics and may present a contamination risk (such as Creutzfeldt-Jakob disease).

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician or practitioner.

Storage and Transport: The valve systems in Integra packaging may be stored at room temperature. They will withstand the conditions of normal transport. *This product is sterilized with ethylene oxide. Do not use if the package is open or damaged. Use the device prior to the «Use by» date on the package label.*

Caution - Do not resterilize. Integra will not be liable for any or all damages including, but not limited to, direct, indirect, incidental, consequential or punitive damages resulting from or related to resterilization.

OSV II PRESSURE / FLOW TEST

Each OSV II Valve System is tested at the end of the manufacturing process. Integra does not recommend testing the valve in the operating room before implantation due to risk of possible contamination and because the conditions required for this testing are difficult to meet in a clinical environment.

PREIMPLANT TESTING FOR PATENCY

Should verification of valve patency be desired prior to implantation, place the inlet of the valve system in a container of sterile apyrogenic saline solution. Connect a syringe, with the Luer connector attached to the drainage catheter, and gently draw the solution through the system to displace air (Figure 3). Remove syringe and Luer connector. Holding the valve's inlet end in the sterile solution, allow the drainage catheter to hang freely (Figure 3). The valve is patent if solution drips from the drainage catheter. Time between each drop will be approximately 4-6 seconds.

Note: This procedure only tests the valve patency. It does not allow verification of the 3 different stages of operation of the Integra OSV II Valve System.

RECOMMENDED IMPLANTATION PROCEDURE

Before implantation, the valve serial number should be recorded in the patient's file.

Note: After implantation, complete the patient identification card for the patient or family.

Clinical Judgment

Incision size, size and type of nonabsorbable suture material, etc. should be selected according to the surgeon's experience and preference. Procedure may be varied in accordance with the surgeon's clinical judgment.

Ventriculoperitoneal Approach

The described procedure is based on the use of a separate ventricular catheter and a specific type of peritoneal introducer. If a valve with an integral ventricular catheter is used, follow the procedure «Introduction of Integral Ventricular Catheters».

1. Place the patient in the supine position. Administer anesthesia. To facilitate tunneling, the head is turned to the contralateral side, placing the abdomen, thorax, neck and skull of the occipital region on a straight plane. Padding may be placed under the neck to facilitate alignment.

Note: The skull incision location depends on the ventricular approach (frontal, temporal or occipital).

2. Expose the skull by making a 25 mm linear incision. Using a perforator, make a burr hole 3.0 mm or larger. Using an electrocoagulator, make an opening in the dura just large enough for ventricular catheter passage. For systems with burr hole cap: using a perforator, make a burr hole larger than 6.5 mm. The technique may be varied according to the surgeon's experience.
3. Using blunt dissection, create a small pocket under the skin to avoid unnecessary pull strength on the drainage catheter when positioning the shunt. Pass tunneler from the subcutaneous pocket to the abdominal area, directing the tunneler tip by trapping it between the thumb and the index finger.

Make the abdominal incision. Pass the catheter from the head to the abdominal incision. Pull on the distal end of the drainage catheter. Position the valve under the skin.

4. Place the introducing rod into the ventricular catheter. If the Right Angle Guide (RAG) is used, introduce the ventricular catheter into the center of the RAG and insert the catheter into the ventricle. Holding the catheter with covered forceps, remove the introducing rod. Verify adequate CSF flow.

For systems with burr hole cap: trim the ventricular catheter to the desired length. Insert the ventricular catheter into the burr hole reservoir and tie with non absorbable suture. Using the introducing rod, insert the catheter into the ventricle. Verify adequate flow.

Connect a syringe with the provided Luer connector to the ventricular catheter and aspirate 2 to 3 ml of CSF to eliminate possible debris from the CSF.

Note: If desired, sample CSF for analysis or inject antibiotics at this time.

- Place the RAG into the burr hole. The channel should point in the direction of the valve. Bend the exposed catheter into the channel of the RAG until the catheter snaps in place.

Note: The use of the RAG in premature infants or neonates is at the discretion of the surgeon. If the RAG is not used, the surgeon should trim the rim of the burr hole where the catheter emerges to provide a smooth angular transition for the catheter.

Trim the ventricular catheter to length (approximately 2 cm from the burr hole).

- Fill the valve system with a sterile apyrogenic saline solution by gentle aspiration through the distal catheter, or using the provided tubing and luer connector for models without distal catheters. Carefully insert the integral inlet straight connector of the valve into the ventricular catheter. Do not dislodge the catheter. Tie the catheter onto the connector and to the underlying fascia, using appropriate sutures.

Note: For configurations with burr hole cap, attach the cap to the reservoir.

- Connect a syringe with the provided Luer connector to the drainage catheter and rinse the system by gently aspirating 2 to 3 ml of CSF to eliminate possible debris and purge air.

- If desired, trim the drainage catheter.

- Introduce the drainage catheter in the peritoneal cavity using a trocar or standard laparotomy technique consistent with the surgeon's experience.

- Caution:** Prior to the introduction of the drainage catheter into the abdominal cavity, it is advisable to confirm the patency of the system by observing the formation of CSF droplets at the end of the drainage catheter.

Notes:

- The CSF droplets formation rate may be lower than that observed with conventional DP shunts because of the flow restriction mode of the OSV II Valve System.

- In children, insert sufficient drainage catheter length into the abdominal cavity to allow for growth.

- Close incisions.

Note: X-ray the complete system just after implantation for future reference to determine whether system components have shifted.

Introduction of Integral Ventricular Catheter (Alternative for Steps 4-5)

If the RAG is used, slide the catheter through the center hole of the RAG. Position the RAG as far along the catheter as possible. Lock the stylet of the ventricular catheter introducer into its cannula as shown in Figure 4. Insert the tip of the introducer into the hole nearest the tip of the ventricular catheter. Align the catheter with the introducer cannula and snap it into the adjustable plastic collar.

Caution: To avoid improper placement, do not overstretch the ventricular catheter when positioning it on the introducer. Puncture the ventricle with this assembly. Carefully withdraw the stylet from the introducer cannula, causing the ventricular catheter tip to disengage from the stylet.

Note: Sample CSF at this time if desired by withdrawing fluid from the cannula.

Disengage the catheter from the collar on the cannula, ensuring that the catheter's position is maintained. Remove the introducer cannula from the burr hole. For instructions on the use of the RAG, refer to Steps 4 and 5 of the «Ventriculoperitoneal Approach» procedure.

Ventriculoatrial Approach

Caution: The procedure is based on the use of a separate ventricular catheter. If a valve with an integral ventricular catheter is used, follow the procedure «Introduction of Integral Ventricular Catheter».

- Position the patient and open the skull as described in Steps 1 and 2 of the «Ventriculoperitoneal Approach» procedure above. Using blunt dissection, create a small pocket under the skin to avoid unnecessary pull strength on the drainage catheter when positioning the shunt.

- Make a neck incision over the point where the angle of the mandible crosses the anterior edge of the sternocleidomastoid muscle. If the internal jugular vein is used, facilitate entry by mobilizing the common facial vein at its junction with the internal jugular vein.

- Tunnel the drainage catheter from the skull subcutaneous pocket to the neck incision. Position the valve under the skin.

- Insert the ventricular catheter following the Steps 4, 5, and 6 of the «Ventriculoperitoneal Approach» procedure.

- Trim the drainage catheter at the neck incision. Fill the distal part with sterile apyrogenic saline solution. Clamp the proximal end. Introduce the drainage catheter into the vein. Position the distal end of the drainage catheter in the heart's right atrium at the level of the 6th or 7th thoracic vertebra. Determine the exact location of the catheter by X-ray, radioscopy, ECG, pressure analysis or echography.

Note: To enhance X-ray visualization, the atrial catheter may be filled with contrast.

- Carefully maintaining the catheter in place, trim its proximal part, and the valve outlet tubing to length. Connect the two together using a straight connector. Tie securely using appropriate sutures.

- Close incisions.

Note - X-ray the complete system just after implantation for future reference to determine whether system components have shifted.

Valve Pumping after Implantation

Percutaneous depression of the antechamber with the finger will force CSF out of the valve in both the proximal and distal directions. Occlude either the distal or proximal side of the antechamber by finger pressure (according to Figure 5), then depress the antechamber. This will push the fluid in the direction opposite to the occlusion.

Warning : Shunt obstruction may occur in any component of a shunt system and should be diagnosed by clinical symptoms and X-rays. Valve pumping testing may not be adequate to diagnose occlusion of catheters.

Injection into the Valve System

The antechamber can be used to access CSF for injections. Use of 25G or smaller beveled needles is recommended. Care must be taken during insertion and removal of the needle to avoid tearing the silicone material. It is also advised to avoid puncturing several times at the same place. Taking into account the above recommendations, the antechamber is qualified to sustain up to 40 punctures without leaking.

Caution: Do not overpressurize the valve system. Excessive flushing pressures may lead to valve damage.

MRI SAFETY INFORMATION

Non-clinical testing has demonstrated that the implanted OSV II Valve systems components are MR CONDITIONAL. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3.0 T or less
- Maximum spatial field gradient of 3500 gauss/cm (35 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4 W/kg (First Level Controlled Mode)

Do not take non-implanted accessories used to aid in introduction of the shunt device such as stylets and catheter guide into the MR environment. They have not been evaluated for safety in the MR environment.

Under the scan conditions defined above, the device is expected to produce a maximum temperature rise of less than 0.4 °C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the OSV II Valve system extends approximately 10 mm from the device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

The flow characteristics of the shunt are not expected to be affected by MR exposure because the valve mechanism does not have any metallic components.

Note: Caution should be exercised when connecting the valve to components that may contain interfering or non-compatible materials.

HOW SUPPLIED

Integra OSV II and OSV II Low Pro Valve Systems are supplied sterile and non-pyrogenic in a double-wrap packaging.

Disposal

After patient use, the system must be handled as biohazardous material and disposed of in accordance with applicable Federal, State, local, or international environmental requirements following facility protocols.

PRODUCT INFORMATION DISCLOSURE

INTEGRA HAS EXERCISED REASONABLE CARE IN THE CHOICE OF MATERIALS AND MANUFACTURE OF THIS PRODUCT. INTEGRA EXCLUDES ALL WARRANTIES, WHETHER EXPRESSED OR IMPLIED BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. INTEGRA SHALL NOT BE LIABLE FOR ANY INCIDENTAL OR CONSEQUENTIAL LOSS, DAMAGE, OR EXPENSE, DIRECTLY OR INDIRECTLY ARISING FROM USE OF THIS PRODUCT. INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY PERSON TO ASSUME FOR IT ANY OTHER OR ADDITIONAL LIABILITY OR RESPONSIBILITY IN CONNECTION WITH THESE PRODUCTS. INTEGRA INTENDS THAT THIS DEVICE SHOULD BE USED ONLY BY PHYSICIANS WITH EDUCATIONAL AND TRAINING BACKGROUND ENABLING THE PROPER USE OF THE DEVICE.

Returned Goods Policy

Authorization from customer service must be obtained prior to returning the product.

Sterile products must be returned in unopened packages with manufacturer's seals intact to be accepted for replacement or credit unless returned due to a complaint or product defect.

Determination of a product defect will be made by Integra.

Credit will be issued for goods returned prior to 90 days from ship date and may be subject to a restocking charge of up to 20%. This assumes that the product returned is not damaged and can be verified to have not been used or opened.

Product Order Information

All products can be ordered through your Integra NeuroSpecialist or customer service representative or by contacting:

Integra LifeSciences Corporation
1100 Campus Road
Princeton, NJ 08540, U.S.A.
Telephone: 1-800-654-2873
Outside the US: +1-609-275-0500
Fax: +1-609-275-5363

Service and Repair
For service, repair or replacements, contact your local Integra representative.

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