



# Patient I.D. Card

**Codman<sup>®</sup> Hakim<sup>®</sup>**

Programmable Valve

In case of medical emergency  
or medical treatment,  
please show this **Patient I.D. Card !**

# Patient Data

## PATIENT DATA

First Name

Last Name

Date of Birth

City

Street

Country

Phone

## EMERGENCY CONTACT

First Name

Last Name

Phone

# Implantation Information

## IMPLANT

Implant date

Implanting hospital

Hospital phone

Doctor name

Doctor phone

Medical condition

## VALVE PLACEMENT

Frontal

Occipital

# Product Information

## PROGRAMMABLE VALVE

The Codman® Hakim® Programmable Valve features 18 pressure settings.

Product Code



Serial Number

- |                                          |                                            |                                |
|------------------------------------------|--------------------------------------------|--------------------------------|
| <input type="checkbox"/> Cylindrical     | <input type="checkbox"/> Regular           | <input type="checkbox"/> Micro |
| <input type="checkbox"/> Micro Burr Hole | <input type="checkbox"/> Regular Burr Hole |                                |

## WARNING

- Use only the Codman® Hakim® Programmer or VPV® System to adjust and indicate the setting of the Codman® Hakim® programmable valve.
- The use of MRI systems (up to 3 Tesla) will not damage the valve mechanism, but may change the operating pressure of the valve. Following any MRI procedure, the clinician should confirm the valve pressure setting and reprogram the valve if necessary.
- The valve setting is adjusted with the application and manipulation of strong magnets. A change to the valve setting is unlikely to occur under normal circumstances. However, magnetic fields should not be placed in close proximity to the valve due to the possibility of an unintentional setting change.
- If these or other symptoms occur contact your Doctor immediately: headaches, vomiting, drowsiness, weakness, irritability or altered mental state.

# Product Information



- 1 Turn lateral radiograph that label 1 is in 1 o'clock position
- 2 Read position of the label according the schema (example shows adjustment 110)

# Post Procedural Data

## POST PROCEDURAL DATA

Date

Pressure Setting

Doctor Notes

## POST PROCEDURAL DATA

Date

Pressure Setting

Doctor Notes

## POST PROCEDURAL DATA

Date

Pressure Setting

Doctor Notes

# Post Procedural Data

## POST PROCEDURAL DATA

Date

Pressure Setting

Doctor Notes

## POST PROCEDURAL DATA

Date

Pressure Setting

Doctor Notes

## POST PROCEDURAL DATA

Date

Pressure Setting

Doctor Notes

## Important

The holder of this Patient Identification Card has a hydrocephalus shunt system implanted for the drainage of cerebrospinal fluid.

In case of lethargy, sleepiness, confusion, fainting, vertigo, epileptic seizure, neck stiffness, fever of unknown origin, damage to the operative sutures, bruising or skin injuries around the operative site, or exposed shunt components, **contact your doctor or hospital immediately.**





## Indications

The Codman® Hakim® Programmable Valves are implantable devices that provide constant intraventricular pressure and drainage of CSF for the management of hydrocephalus.

## Contraindications

The Codman® Hakim® Programmable Unitized Valve Systems are not recommended for atrial placement. Use the non-unitized versions for this procedure. These devices are contraindicated in patients receiving anticoagulants or known to have a bleeding diathesis. Avoid shunt implantation if infection is present within the body. Delay the shunt procedure when infections such as meningitis, ventriculitis, peritonitis, bacteremia, and septicemia are present.

---

### For more information or to place an order, please contact:

Sales & Marketing EMEA  
Integra LifeSciences Services (France) SAS  
Immeuble Séquoia 2 97 allée A. Borodine  
Parc technologique de la Porte des Alpes  
69800 Saint Priest - FRANCE  
Phone: +33 (0)4 37 47 59 00  
Fax: +33 (0)4 37 47 59 99



Integra LifeSciences Switzerland Sarl  
Rue Girardet 29 (2nd Floor)  
Le Locle ■ Neuchâtel CH-2400  
Switzerland

Codman, Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. Hakim is a registered trademark of Hakim USA, LLC and is used under license.

©2020 Integra LifeSciences Corporation. All rights reserved. Last revision date: 05/2020.  
1077914-2-EN