


Codman[®]
SPECIALTY SURGICAL

A DIVISION OF INTEGRA LIFESCIENCES



Integra[®] **Codman**
Neuro Accessories
Product Catalog

Customer Service and warranty information	3
Access & hemostasis	4
Scalp clips	4
Scalp clips appliers	4
Disposable perforators	5
Craniotomy kit	5
Protect the brain	6
Surgical patties	6
Surgical strips	7
Bicol Collagen sponge	7
Cranioplastic Kit	8
Indications & Contraindications	9
Products References	10
Actual size drawings	12
CODMAN Surgical Patties sizes	12
CODMAN Surgical Strips sizes	13

Customer Services:

International:

+33 (0)437 47 59 50
 +33 (0)437 47 59 25 (Fax)
 csemea@integralife.com

Germany:

+49 (0)2102 5535 6200
 +49 (0)2102 5536 636 (Fax)
 custsvcgermany@integralife.com

Switzerland:

+41 (0)2 27 21 23 00
 +41 (0)2 27 21 23 99 (Fax)
 custsvcsuisse@integralife.com

Austria:

+43 (0)720816067
 +43 (0)19287201 (Fax)
 custsvcaustria@integralife.com

Ireland:

+353 1800 901567
 +353 1822 5952 (Fax)
 custsvcire@integralife.com

United Kingdom:

+44 (0)1264 312 725
 +44 (0)1264 312 821 (Fax)
 custsvcs.uk@integralife.com

Belgium & Luxembourg:

+32 (0)2 257 4130
 +32 (0)2 253 2466 (Fax)
 custsvcbenlux@integralife.com

Italy:

+39 (0)2 577 89 21
 +39 (0)2 575 113 71 (Fax)
 custsvcitaly@integralife.com

France :

+33 (0) 437 47 59 10
 +33 (0) 437 47 59 29 (Fax)
 custservfrance@integralife.com

Netherlands:

+31 (0)852083167
 +31 (0)207093627 (Fax)
 custsvcnetherlands@integralife.com

Integra Limited Warranty

INTEGRA LIFESCIENCES CORPORATION and its wholly owned subsidiaries («INTEGRA») warrant to INTEGRA authorized distributors and the original purchaser only that each new INTEGRA or CODMAN product is free from manufacturing defects in material and workmanship under normal use and service from the date of delivery by INTEGRA (or its authorized distributor) to the original purchaser, but in no event beyond the expiration date stated on any product labeling. For purposes of products sold by INTEGRA through an authorized distributor of INTEGRA, «original purchaser» shall include the purchaser of INTEGRA products to whom the distributor first sells the product.

- Surgical instruments are guaranteed to be free from defects in material and workmanship when maintained and cleaned properly and used normally for their intended purpose.
- Any covered product that is placed by INTEGRA under a lease, rental or installment purchase agreement and that requires repair service during the term of such placement agreement shall be repaired in accordance with the terms of such agreement.

If any covered defect occurs during the warranty period or term of such placement agreement, the purchaser or distributor should communicate directly with INTEGRA. If purchaser or distributor seeks to invoke the terms of this warranty, the product must be returned to INTEGRA. The defective product should be returned promptly, properly packaged and postage prepaid. Loss or damage in return shipment to INTEGRA shall be at sender's risk. INTEGRA's sole responsibility under this warranty shall be repair or replacement, at INTEGRA's sole discretion at INTEGRA's expense, subject to the terms of this warranty and applicable agreements.

IN NO EVENT SHALL INTEGRA BE LIABLE FOR ANY INCIDENTAL, INDIRECT, CONSEQUENTIAL OR PUNITIVE DAMAGES IN CONNECTION WITH THE ACQUISITION OR USE OF ANY INTEGRA PRODUCT. Further, this warranty shall not apply to, and INTEGRA shall not be responsible for, any loss arising in connection with the purchase or use of any INTEGRA product that has been repaired by anyone other than an authorized INTEGRA service representative or altered in any way so as, in INTEGRA's judgment, to affect its stability or reliability, or which has been subject to misuse, negligence or accident, or which has been used otherwise than in accordance with the instructions furnished by INTEGRA.

THIS INTEGRA LIMITED WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON INTEGRA'S PART OR THE PART OF ITS DISTRIBUTORS, AND INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY REPRESENTATIVE OR OTHER PERSON TO ASSUME FOR IT ANY OTHER LIABILITY IN CONNECTION WITH INTEGRA'S PRODUCTS.

INTEGRA DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION OR WARRANTY OF QUALITY AS WELL AS ANY EXPRESS OR IMPLIED WARRANTY TO PATIENTS. No warranty or guarantee may be created by any act or statement nor may this Standard Warranty be modified in any way, except as a result of a writing signed by an officer of INTEGRA. These limitations on the creation or modification of this warranty may not be waived or modified orally or by any conduct.

IN NO EVENT SHALL INTEGRA AUTHORIZED DISTRIBUTORS BE LIABLE TOWARDS THE ORIGINAL PURCHASER FOR ANY INCIDENTAL, INDIRECT, CONSEQUENTIAL OR PUNITIVE DAMAGES IN CONNECTION WITH THE ACQUISITION OR USE OF ANY INTEGRA PRODUCT. Further, this warranty shall not apply to, and INTEGRA authorized distributors shall not be responsible towards the original purchaser for, any loss, arising in connection with the purchase or use of any INTEGRA product that has been repaired by anyone other than an authorized INTEGRA service representative or altered in any way so as to affect its stability or reliability, or which has been subject to misuse, negligence or accident, or which has been used otherwise than in accordance with the instructions furnished by INTEGRA. THIS INTEGRA DISTRIBUTOR LIMITED WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES TOWARDS THE ORIGINAL PURCHASER, EXPRESS OR IMPLIED, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES TOWARDS THE ORIGINAL PURCHASER ON INTEGRA AUTHORIZED DISTRIBUTOR'S PART.

INTEGRA AUTHORISED DISTRIBUTORS DISCLAIM ALL OTHER WARRANTIES TOWARDS THE ORIGINAL PURCHASER, EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION OR WARRANTY OF QUALITY AS WELL AS ANY EXPRESS OR IMPLIED WARRANTY TO PATIENTS.

Scalp clips

Scalp clips to provide hemostasis to a scalp wound edge

The waved atraumatic design and the adjusted closure tension prevent skin necrosis. A low profile design limits the volume around the operating site. Winged extremities allow an easier handling with the applier.



LEROY Scalp Clips

Reference	Size	Packaging	Compatibility	Additional features
20-1032	Adult	10 clips/package 20 packages/case	Sterile 20-5117 LEROY Clip applying forcep.	Disposable, plastic
20-1035	Infant	10 clips/package 20 packages/case	Sterile 20-5120CN LEROY Clip applying Forcep.	Disposable, plastic



RANEY Scalp clips

Reference	Packaging	Compatibility	Additional features
20-1025	12 clips/package non sterile	20-5115CN RANEY clip applying forcep.	Reusable, stainless steel
20-1037	10 clips/package 20 packages/case	Sterile 20-5115CN RANEY clip applying forcep	Disposable, plastic



LEROY RANEY Scalp clips

Reference	Packaging	Compatibility	Additional features
20-1034	10 clips/package 20 packages/case	Sterile 20-5117 LEROY clip applying forcep.	Disposable, plastic

Scalp clip appliers

Scalp clip appliers and forceps allow for easy application and removal of clips. Both disposable appliers and reusable scalp clip applier forceps are available.

CODMAN DISPOSABLE SCALP CLIP APPLIER, sterile, single use system designed to permit rapid application of scalp clips to the scalp flap for hemostasis during a craniotomy procedure. The disposable applier is designed for one hand use with 12 pre-loaded Leroy-Raney design clips which are easily removed with specially designed removal forceps.

CODMAN Disposable Scalp clip removal and applier forceps

Reference	Packaging	Description
20-1038	10 per box	CODMAN Disposable Scalp clip applier
20-1039		CODMAN Scalp clip removal forceps for 20-1038 clips



LEROY Reusable Scalp clip applying forceps

Reference	Size	Compatibility
20-5117	Adult - 7 1/4 (184mm)	20-1032 or 20-1034 LEROY clips
20-5120CN	Infant - 7 1/4 (184mm)	20-1035 Infant LEROY clips



RANEY Reusable Scalp clip applying forceps

Reference	Size	Compatibility
20-5115CN	6 1/4 (159mm)	20-1037 or 20-1025 RANEY clips



Disposable perforators

Sterile disposable perforators designed for use in perforating the cranium

Safety disengagement mechanism designed to declutch automatically and prevent damage to the dura mater. Supplied sterile and single use, each one tested individually for optimal performance and safety. Available in 3 sizes with colour code to help identify them during surgery.

Reference	Size	Description
26-1221	14mm/11mm	Disposable Perforator
26-1222	11mm/8mm	Disposable Perforator
26-1223	9mm/6mm	Disposable Perforator
26-1246		Cranio blade
26-1247		Wire Pass drill



Craniotomy kit

The craniotomy kit contains sterile and disposable items.

Reference	Quantity	Description
26-1230	1 unit	CODMAN perforators
26-1230	1 unit	Cranio blade
26-1230	1 unit	Wire pass drill



Surgical patties and strips

Codman surgical patties and surgical strips protect brain tissues and other tissues of the central nervous system to meet the demand of delicate intracranial surgery.

Protection

Codman Patties and Strips absorb more than five times their weight in less than a second to help maintain a moist and clean operative site. The Codman Cottonoid® material is highly specialised to provide a soft, porous and conformable protection for delicate tissue. It allows a safe suction through the patties and strips.

Reliability


Cottonoid® non-woven fabric resists to shedding, stretching or linting while providing a cotton-like softness. No debris will remain in the operating site preventing arachnoiditis or peridural fibrosis.

Safety

Tantalum marking makes patties and strips detectable on X-Rays. Sutureweld®'s strong bond holds the indicator string to the patty.



Surgical patties



Real size drawings on page 22.


CODMAN Surgical Patties, sterile, X-ray detectable

Reference	Size	Packaging	Additional features
80-1396	5/16" (7mm)	10 patties per package 20 packages per dispenser box	sterile
80-1397	1/4" x 1 1/2" (7mm x 38mm)	10 patties per package 20 packages per dispenser box	sterile
80-1398	1/4" x 3" (7mm x 76mm)	10 patties per package 20 packages per dispenser box	sterile
80-1399	1/4" x 1/4" (7mm x 7mm)	10 patties per package 20 packages per dispenser box	sterile
80-1400	1/2" x 1/2" (13mm x 13mm)	10 patties per package 20 packages per dispenser box	sterile
80-1401	3/4" x 3/4" (19mm x 19mm)	10 patties per package 20 packages per dispenser box	sterile
80-1402	1/2" x 1" (13mm x 25mm)	10 patties per package 20 packages per dispenser box	sterile
80-1403	1" x 1" (25mm x 25mm)	10 patties per package 20 packages per dispenser box	sterile
80-1404	1/2" x 1 1/2" (13mm x 38mm)	10 patties per package 20 packages per dispenser box	sterile
80-1406	1/2" x 2" (13mm x 52mm)	10 patties per package 20 packages per dispenser box	sterile
80-1407	1/2" x 3" (13mm x 76mm)	10 patties per package 20 packages per dispenser box	sterile
80-1408	1" x 3" (25mm x 76mm)	10 patties per package 20 packages per dispenser box	sterile
80-1409	3" x 3" (76mm x 76mm)	10 patties per package 20 packages per dispenser box	sterile

Surgical strips

CODMAN Surgical Strips, X-ray detectable

Reference	Size	Packaging	Additional features
80-1449	1/8" x 6" (3 x 152mm)	10 strips per package 20 packages per dispenser box	sterile
80-1450	1/4" x 6" (6 x 152mm)	10 strips per package 20 packages per dispenser box	sterile
80-1451	1/2" x 6" (13 x 152mm)	10 strips per package 20 packages per dispenser box	sterile
80-1452	3/4" x 6" (19 x 152mm)	10 strips per package 20 packages per dispenser box	sterile
80-1453	1" x 6" (25 x 152mm)	10 strips per package 20 packages per dispenser box	sterile
80-1454	1 1/2" x 6" (38 x 152mm)	10 strips per package 20 packages per dispenser box	sterile
80-1455	2" x 6" (51 x 152mm)	10 strips per package 20 packages per dispenser box	sterile
80-1456	3" x 6" (76 x 152mm)	10 strips per package 20 packages per dispenser box	sterile
80-1457	3 1/2" x 6" (89 x 152mm)	10 strips per package 20 packages per dispenser box	sterile



Real size drawings on page 23.

Bicol® Collagen sponge

The BICOL Collagen Sponge is a moistening agent and protective device. One surface is smooth so that once moistened it slides over the brain conforming to the area to be protected. The opposite surface has sufficient grain to provide friction to a retractor blade or surgical pattie for manipulation. The BICOL sponge buffers the pressure of retractor blades.

BICOL Sponges

Reference	Size	Packaging	Additional features
80-1460	3"x4" (76 x 102 mm)	Individually, sold 10 packages per box	sterile



CRANIOPLASTIC™ Kit

CRANIOPLASTIC™ Kit

Reference	Packaging
42-1280	Two complete sets, consists of two bags sterile powder (net weight 30 grams each), two vials sterile liquid (net contents 17 milliliters each) and two sterile matrix bags.



Bicol® Collagen sponge

Indications: The Bicol® collagen sponge is for use in neurological procedures requiring prolonged retraction and exposure of the brain

Contraindications: The Bicol® collagen sponge should not be packed into vascular sheaths, neural tracts, interfascial planes, sinuses, foramina, or areas of confined anatomical structure. If confined packing should become necessary in certain situations, care should be taken to avoid the possible development of undue pressure on vital structures.

Cranioplastic

Indications: CRANIOPLASTIC is a resinous material for repairing cranial defects.

Contraindications: CRANIOPLASTIC resin is contraindicated in cases of active infection. CRANIOPLASTIC resin should not be used in patients previously sensitized to methyl methacrylate.

Patties & Strips

Indications: The surgical patties and surgical strips are indicated for the use in protection of tissue, including brain and other tissues of the central nervous system, during surgery.

Contraindications: Do not use for other indications than the ones specified in the instructions for use.

Perforators

Indications: The CODMAN disposable perforator is for use in perforating the cranium. When properly used, it is designed to automatically disengage once perforation is accomplished and when pressure is removed from the drill point.

Contraindications: Do not use for other indications than the ones specified in the instructions for use.

Raney scalp clips

Indications: The Raney scalp clips provide hemostasis to a scalp wound edge. The Raney scalp clip applying forceps are recommended for applying the Raney scalp clips.

Contraindications: Do not use for other indications than the ones specified in the instructions for use.

Numerical Index

20-1025	Raney scalp clips - reusable	p.4
20-1032	Leroy scalp clips - Adult.....	p.4
20-1034	Leroy-Raney scalp clips	p.4
20-1035	Leroy scalp clips - Infant.....	p.4
20-1037	Raney scalp clips - Disposable.....	p.4
20-1038	Disposable scalp clip applier	p4
20-1039	Disposable scalp clip removal forceps.....	p.4
20-5115CN	Raney scalp clip applying forceps	p.4
20-5117	Leroy scalp clip applying forceps - Adult	p.4
20-5120CN	Leroy scalp clip applying forceps - Infant	p.4
26-1221	14mm disposable perforator	p.5
26-1222	11mm disposable perforator.....	p.5
26-1223	9mm disposable perforator	p.5
26-1230	Craniotomy kit	p.5
26-1246	Cranio blade	p.5
26-1247	Drill, wire pass.....	p.5
42-1280	Cranioplastic™.....	p.8
80-1396	Micro patties round 5/16" (.79 cm).....	p.6, p.12
80-1397	Patties 1/4" x 1 1/2" (64cm x 3.81cm).....	p.6, p.12
80-1398	Patties 1/4" x 3" (64cm x 7.62cm)	p.6, p.12
80-1399	Patties 1/4" x 1/4" (64cm x 64cm)	p.6, p.12
80-1400	Patties 1/2" x 1/2" (1.27cm x 1.27cm).....	p.6, p.12
80-1401	Patties 3/4" x 3/4" (1.91cm x 1.91cm).....	p.6, p.12
80-1402	Patties 1/2" x 1" (1.27cm x 2.54cm)	p.6, p.12
80-1403	Patties 1" x 1" (2.54cm x 2.54cm).....	p.6, p.12
80-1404	Patties 1/2" x 1 1/2" (1.27cm x 3.81cm)	p.6, p.12
80-1406	Patties 1/2" x 2" (1.27cm x 5.08cm)	p.6, p.12
80-1407	Patties 1/2" x 3" (1.27cm x 7.62cm)	p.6, p.12
80-1408	Patties 1" x 3" (2.54cm x 7.62cm)	p.6, p.12
80-1409	Patties 3" x 3" (7.62cm x 7.62cm)	p.6, p.12
80-1449	Strips 1/8" x 6" (0.3 cm x 15.2 cm).....	p.7, p.13
80-1450	Strips 1/4" x 6" (0.6 cm x 15.2 cm)	p.7, p.13
80-1451	Strips 1/2" x 6" (1.3 cm x 15.2 cm).....	p.7, p.13
80-1452	Strips 3/4" x 6" (1.9 cm x 15.2 cm)	p.7, p.13
80-1453	Strips 1" x 6" (2.5 cm x 15.2 cm)	p.7, p.13
80-1454	Strips 1 1/2" x 6" (3.8cm x 15.2cm)	p.7, p.13
80-1455	Strips 2" x 6" (5.08 cm x 15.2cm)	p.7, p.13
80-1456	Strips 3" x 6" (7.6cm x 15.2cm)	p.7, p.13
80-1457	Strips 3 1/2" x 6" (8.9cm x 15.2cm)	p.7, p.13
80-1460	BICOL Collagen Sponge	p.7

Bicol[®] collagen sponge

80-1460	BICOL Collagen Sponge	P.7
---------	-----------------------------	-----

CRANIOPLASTIC[™] Kit

42-1280	CRANIOPLASTIC [™] Kit	P.8
---------	--------------------------------------	-----

Patties and strips – Packed per 200 per box

80-1396	Micro patties round 5/16" (.79 cm)	p.6, p.12
80-1397	Patties 1/4" x 11/2" (64cm x 3.81cm)	p.6, p.12
80-1398	Patties 1/4" x 3" (64cm x 7.62cm)	p.6, p.12
80-1399	Patties 1/4" x 1/4" (64cm x 64cm)	p.6, p.12
80-1400	Patties 1/2" x 1/2" (1.27cm x 1.27cm)	p.6, p.12
80-1401	Patties 3/4" x 3/4" (1.91cm x 1.91cm)	p.6, p.12
80-1402	Patties 1/2" x 1" (1.27cm x 2.54cm)	p.6, p.12
80-1403	Patties 1" x 1" (2.54cm x 2.54cm)	p.6, p.12
80-1404	Patties 1/2" x 1 1/2" (1.27cm x 3.81cm)	p.6, p.12
80-1406	Patties 1/2" x 2" (1.27cm x 5.08cm)	p.6, p.12
80-1407	Patties 1/2" x 3" (1.27cm x 7.62cm)	p.6, p.12
80-1408	Patties 1" x 3" (2.54cm x 7.62cm)	p.6, p.12
80-1409	Patties 3" x 3" (7.62cm x 7.62cm)	p.6, p.12
80-1449	Strips 1/8" x 6" (0.3 cm x 15.2 cm)	p.7, p.13
80-1450	Strips 1/4" x 6" (0.6 cm x 15.2 cm)	p.7, p.13
80-1451	Strips 1/2" x 6" (1.3 cm x 15.2 cm)	p.7, p.13
80-1452	Strips 3/4" x 6" (1.9 cm x 15.2 cm)	p.7, p.13
80-1453	Strips 1" x 6" (2.5 cm x 15.2 cm)	p.7, p.13
80-1454	Strips 1 1/2" x 6" (3.8cm x 15.2cm)	p.7, p.13
80-1455	Strips 2" x 6" (5.08 cm x 15.2cm)	p.7, p.13
80-1456	Strips 3" x 6" (7.6cm x 15.2cm)	p.7, p.13
80-1457	Strips 3 1/2" x 6" (8.9cm x 15.2cm)	p.7, p.13

Perforators

26-1221	14mm disposable perforator	p.5
26-1222	11mm disposable perforator	p.5
26-1223	9mm disposable perforator	p.5
26-1230	Craniotomy kit	p.5
26-1246	Cranio blade	p.5
26-1247	Drill, wire pass	p.5


Scalp clips

20-1025	Raney scalp clips - reusable	P.4
20-1032	Leroy scalp clips - Adult	P.4
20-1034	Leroy-Raney scalp clips	P.4
20-1035	Leroy scalp clips - Infant	P.4

20-1037	Raney scalp clips - Disposable	P.4
20-1038	Disposable scalp clip applicator	P.4
20-1039	Disposable scalp clip removal forceps	P.4
20-5115CN	Raney scalp clip applying forceps	P.4
20-5117	Leroy scalp clip applying forceps - Adult	P.4
20-5120CN	Leroy scalp clip applying forceps - Infant	P.4


CODMAN Surgical Patties, sterile, X-ray detectable

Reference	Size		Packaging		Additional features
	in inch	in mm	units/pack	packs/box	
80-1396	5.16"	7	10	20	sterile
80-1397	1.4" x 11.2"	7x38			
80-1398	1.4" x 3"	7x76			
80-1399	1.4" x 1.4"	7x76			
80-1400	1.2" x 1.2"	13x13			
80-1401	3.4" x 3.4"	19x19			
80-1402	1.2" x 1"	13x25			
80-1403	1" x 1"	25x25			
80-1404	1.2" x 11.2"	13x38			
80-1406	1.2" x 2"	13x52			
80-1407	1.2" x 3"	13x76			
80-1408	1" x 3"	25x76			
80-1409	3" x 3"	76x76			



CODMAN Surgical Strips, X-ray detectable



Reference	Size		Packaging		Additional features
	in inch	in mm	units/pack	packs/box	
80-1449	1/8" x 6"	3X152	10	20	sterile
80-1450	1/4" x 6"	6X152			
80-1451	1/2" x 6"	13X152			
80-1452	3/4" x 6"	19X152			
80-1453	1" x 6"	25X152			
80-1454	1 1/2" x 6"	38X152			
80-1455	2" x 6"	51X152			
80-1456	3" x 6"	76X152			
80-1457	3 1/2" x 6"	89X152			



Codman®

SPECIALTY SURGICAL

A DIVISION OF INTEGRA LIFESCIENCES



Distributed by

Integra LifeSciences Services (France) SAS
Sales & Marketing EMEA
Immeuble Séquoia 2 • 97 allée Alexandre Borodine
Parc technologique de la Porte des Alpes
69800 Saint Priest • FRANCE
+33 (0)4 37 47 59 00 phone • +33 (0)4 37 47 59 99 fax
integralife.eu


Customer Services

International: +33 (0)437 47 59 50 • +33 (0)437 47 59 25 (Fax) • csmea@integralife.com
Austria: +43 (0)720816067 • +43 (0)19287201 (Fax) • custsvcaustria@integralife.com
Belgium & Luxembourg: +32 (0)2 257 4130 • +32 (0)2 253 2466 (Fax) • custsvcbelux@integralife.com
France: +33 (0) 437 47 59 10 • +33 (0) 437 47 59 29 (Fax) • custservfrance@integralife.com
Germany: +49 (0)2102 5535 6200 • +49 (0)2102 5536 636 (Fax) • custsvgermany@integralife.com
Ireland: +353 1800 901567 • +353 1822 5952 (Fax) • custsvcire@integralife.com
Italy: +39 (0)2 577 89 21 • +39 (0)2 575 113 71 (Fax) • custsvitaly@integralife.com
Netherlands: +31 (0)852083167 • +31 (0)207093627 (Fax) • custsvcnetherlands@integralife.com
Switzerland: +41 (0)2 27 21 23 00 • +41 (0)2 27 21 23 99 (Fax) • custsvcsuisse@integralife.com
United Kingdom: +44 (0)1264 312 725 • +44 (0)1264 312 821 (Fax) • custsvcs.uk@integralife.com


©2020 Integra LifeSciences Corporation. All rights reserved. 116236-5-EN
DOCUMENT FOR USE IN EUROPE, MIDDLE-EAST and AFRICA ONLY. Last revision date: 01/2021.

 Integra LifeSciences Production Corporation
11 Cabot Boulevard
Mansfield • MA 02048 • USA



 Integra LifeSciences Switzerland Sarl
Rue Girardet 29 (2nd Floor)
Le Locle • Neuchâtel CH-2400 • Switzerland



 Integra LifeSciences Services (France)
Immeuble Séquoia 2
97 Allée Alexandre Borodine
Parc Technologique de la Porte des Alpes
69800 Saint Priest • France

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region. All text and diagrams contained in this catalog are intended as guidelines and are for illustrative purposes only.
Information included in this catalog is not intended to replace the Directions for Use packaged with each Integra product. Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instruction for use. Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality. WARNING: Applicable laws restrict these products to sale by or on the order of a physician. Codman are registered trademarks of its owners. Bicol, Cottonoid, Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. Products mentioned in this document are CE class I, IIa and III devices. Please contact Integra customer service should you need any additional information on devices classification. All the references numbers mentioned on this document are CE marked according to European council directive 93/42/EEC on medical devices and its relatives, unless specifically identified as "NOT CE MARKED."