

Integra® Jarit®

Sterilization Containers

Limit uncertainty by protecting
your instrument investment



INTEGRA®
LIMIT UNCERTAINTY

Summary of Integra Sterilization Container Product Features

Gaskets within retention plates and container perimeter help ensure a secure seal of the container, lid and filter.

Snap on filter plate with a locking mechanism for changing disposable filters

Round disposable filters with centering hole to assist in proper alignment.

Lid and Bottom manufactured using anodized aluminum alloy. Easy to maintain, corrosion resistant, lightweight yet durable, and conducts heat for material drying period.

4 ID Tag Windows (2 front & 2 back)

ID tags available in 6 colors and can be engraved with up to 18 characters per line.


Low profile latch securely closes lid over container bottom.

Tamper Evident Lock – helps identify sterility breach.

Fold-down handles extend and stop at 90° to aid lifting.

Corner design allows for same size container stacking in storage.

2 Indicator Card Windows (1 front & 1 back). Indicator card has visible bar for clear indication that steam sterilization has been achieved.



One-Half Size (11¼" x 11¾") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4⅛"	730-400	730-400B	730-410	730-410B
5¼"	730-401	730-401B	730-411	730-411B
5⅞"	730-402	730-402B	730-412	730-412B
7¾"	730-403	730-403B	730-413	730-413B
10¼"	730-404	730-404B	730-414	730-414B

Note: Complete container assembly only available in silver.
See below for other color options with lids and bottoms sold separately.



One-Half Size Colored Lids (bottoms sold separately)

Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730-420
Red	730-422
Gold	730-423
Green	730-424
Black	730-425



One-Half Size Wire or Perforated Baskets

Height	Wire Basket	Perforated Basket
2"	740-400	750-400
2¾"	740-401	750-401
4"	740-402	750-402
Lid	-	750-420



Wire Basket



Perforated Basket

Note: The wire baskets come with 'screw-on' feet and the perforated baskets come with 'snap-on' feet.

Three-Quarter Size (18¼" x 11¼") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4¼"	730-500	730-500B	730-510	730-510B
5¼"	730-501	730-501B	730-511	730-511B
5⅞"	730-502	730-502B	730-512	730-512B



Three-Quarter Size Colored Lids (bottoms sold separately)

Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730-520
Red	730-522
Gold	730-523
Green	730-524
Black	730-525




Three-Quarter Size Wire or Perforated Baskets

Height	Wire Basket	Perforated Basket
2"	740-500	750-500
2¾"	740-501	750-501
4"	740-502	750-502
Lid	-	750-520



Wire Basket



Perforated Basket

Note: Baskets are constructed with solid wire mesh to allow ultrasonic cleaning; handles can be positioned to assist with a secure grasp for removal under sterile conditions.

Full-Size (23¼" x 11¼") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4⅛"	730-600	730-600B	730-610	730-610B
5¼"	730-601	730-601B	730-611	730-611B
5⅞"	730-602	730-602B	730-612	730-612B
7¾"	730-603	730-603B	730-613	730-613B
10¼"	730-604	730-604B	730-614	730-614B

Note: Complete container assembly only available in silver.
See below for other color options with lids and bottoms sold separately.



Full-Size Colored Lids (bottoms sold separately)

Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730-620
Red	730-622
Gold	730-623
Green	730-624
Black	730-625




Full-Size Wire or Perforated Baskets

Height	Wire Basket	Perforated Basket
2"	740-600	750-600
2¾"	740-601	750-601
4"	740-602	750-602
Lid	-	750-620



Wire Basket



Perforated Basket

Wide-Body (23½" x 15¾") Containers

- Vacuum Sterilizer Containers - Solid/Perforated Bottom
- Gravity Sterilizer Containers - Perforated Bottom

Height	Solid Bottom Full Container	Solid Bottom Only	Perforated Bottom Full Container	Perforated Bottom Only
4¾"	730-701	730-701B	730-711	730-711B
7"	730-702	730-702B	730-712	730-712B



Wide-Body Lid (bottoms sold separately)

Includes Filter Retention Plate and Latches

Color	Lid Only
Silver	730-720




Wide-Body Wire Baskets

Height	Wire Basket
2"	740-700
4"	740-701



Instrument Stringers and Holders

Catalog Number	Description	Size
117-240	Weinstein Instrument Holder	6"
117-242	Weinstein Instrument Holder	8"
117-244	Weinstein Instrument Holder	10"
117-246	Weinstein Instrument Holder	12"
117-247	Weinstein Instrument Holder	14"
117-248	Weinstein Instrument Holder	16"
117-249	Weinstein Instrument Holder	18"
117-250	Jarit Instrument Holder	6"
117-252	Jarit Instrument Holder	8"
117-254	Jarit Instrument Holder	10"
117-256	Jarit Instrument Holder	12"
117-258	Jarit Instrument Holder	14"
117-260	Jarit Instrument Holder	16"
117-262	Jarit Instrument Holder	18"
117-270	Instrument Stringer	8"
117-272	Instrument Stringer	10"
117-274	Instrument Stringer	12"
117-276	Instrument Stringer	14"
117-278	Instrument Stringer	16"
117-280	Instrument Stringer	18"



Weinstein Instrument Holder

Jarit Instrument Holder

Instrument Stringer

Note: For more instruments stringers please refer to the Jarit® catalog.

Identification Tags

Each tag allows for 2 lines of custom text, 18 characters per line.


Color	Catalog Number
Silver	775-100
Blue	775-101
Red	775-102
Gold	775-103
Green	775-104
Black	775-105



Snap-on Basket Tags (with key ring)

Each small tag allows for 3 lines of custom text, 32 characters per line.
 Each large tag allows for 5 lines of custom text, 32 characters per line.

Color	Small	Size	Large	Size
Silver	775-120	2 3/8" x 1/2"	775-130	2 3/8" x 1 1/4"
Gold	775-121	2 3/8" x 1/2"	775-131	2 3/8" x 1 1/4"
Red	775-122	2 3/8" x 1/2"	775-132	2 3/8" x 1 1/4"
Blue	775-123	2 3/8" x 1/2"	775-133	2 3/8" x 1 1/4"
Green	775-124	2 3/8" x 1/2"	775-134	2 3/8" x 1 1/4"
Black	775-125	2 3/8" x 1/2"	775-135	2 3/8" x 1 1/4"




Tamper Evident Locks and Indicator Cards

Catalog Number	Description
775-400	Tamper Evident Locks, single use, blue, box of 1000
775-300	Indicator Cards, box of 250




Tamper Evident Lock




Indicator Card

Retention Plate and Filters

Catalog Number	Description
775-212-INT	Retention Plate for Paper Filters, Round
775-200	Disposable Filters for Regular Containers PKG/1000



Retention Plate




Disposable Filters

Instrument Rack

- 9 slots
- 5mm - 10mm diameter silicone bars
- Compatible with the 8" Full-Size container
- The rack fits inside the 2" or 3" wire or perforated basket, allowing storage underneath for other instruments and auxiliary items

Catalog Number	Description
790-010	Instrument Rack



Silicone Holders

Catalog Number	Description	Size
775-630	Silicone Holder #1	8mm x 26mm x 40mm
775-631	Silicone Holder #2	8mm x 26mm x 40mm
775-632	Silicone Holder #3	8mm x 26mm x 40mm
775-650	Silicone Holder Clamps, 4 slot, plastic	-



#1




#2



#3



Silicone Mats

Catalog Number	Description
776-400	One-Half Size Silicone Mat
776-500	Three-Quarter Size Silicone Mat
776-600	Full-Size Silicone Mat



Partition Sheets and Holding Pins

Catalog Number	Description	Size
775-610	Partition Sheet	2" x 7/8" (50mm x 20mm)
775-611	Partition Sheet	5" x 7/8" (130mm x 20mm)
775-612	Partition Sheet	9 7/8" x 7/8" (225mm x 20mm)
775-613	Partition Sheet	18 1/8" x 7/8" (460mm x 20mm)
775-600	Holding Pin for Wire Baskets	
775-601	Holding Pin for Partition Sheets	



Holding Pin for Wire Baskets



Holding Pin for Partition Sheets

Holding Clamps

Catalog Number	Size Range
775-620	1/8" to 1/2" (7mm-12mm)
775-621	1/2" to 5/8" (12mm-16mm)
775-622	5/8" to 1" (16mm-26mm)
775-623	1" to 1 1/2" (26mm-36mm)
775-624	1 1/2" to 1 3/4" (36mm-45mm)



Screw-on Feet

Catalog Number	Description
780FT-100	Plastic Foot for Baskets, Pkg/1

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

integralife.com/contact

