

Mini-Port Laparoscopic set
Integra® MicroFrance®

3.5 mm Laparoscopic Instruments

For mini-invasive surgery



Integra® MicroFrance®

3.5 mm Laparoscopic Instruments

The 3.5 mm Mini-Port laparoscopic set for mini-invasive surgery


■ Handcrafted in France

MicroFrance, thanks to more than 45 years of experience in manufacturing surgical instruments, is utilizing his expertise to develop a comprehensive offer of mini-invasive laparoscopic instruments.



■ Mini-invasive surgery

With a 3.5 mm diameter, the MicroFrance mini-port set enables a reduced incision size¹ for better patient recovery².



■ Quality without compromise

The diameter reduction still guarantees a high instruments rigidity. The surgical gesture remains precise thanks to the forceps tactile feedback.

■ Versatility & ergonomics

A large offer of both 1-piece and semi-modular instruments, with multiple handle shapes for a flexible prehension.

1- Compared to conventional laparoscopy, with 5 and 10 mm incisions

2- Novitsky YW, Kercher KW, Czerniach DR, et al. Advantages of mini-laparoscopic vs conventional laparoscopic cholecystectomy: results of a prospective randomized trial. Arch Surg. 2005;140(12):1178–1183

Integra® MicroFrance®

3.5 mm Laparoscopic Instruments

Mini-Port Instruments family

A large product portfolio available in 2 lengths: standard and pediatric, with inserts offering rigidity and enabling precise surgical gestures

Babcock forceps



Inserts*:




CEV391IU: 335 mm
CEV891IU: 235 mm

1-piece:



CEV391BU: 355 mm
Montupet CEV891MU: 220 mm

Alligator jaws forceps



Inserts*:



CEV396IU: 330 mm
CEV896IU: 230 mm


1-piece:



De Lagausie CEV997MU: 215 mm

*Inserts compatible with every rotative handles available on page 8 & 9

Fenestrated grasping forceps




Inserts*:



CEV395IU: 335 mm
CEV895IU: 235 mm

1-piece:



CEV395BU: 355 mm
Montupet CEV895MU: 215 mm

1-piece:



De Lagausie CEV995MU: 220 mm

Grip forceps




Inserts*:



CEV399IU : 335 mm
CEV899IU: 220 mm

1-piece:




CEV399BU: 355 mm
CEV899U: 220 mm

Integra® MicroFrance®

3.5 mm Laparoscopic Instruments

Mini-Ports Instruments family

Dissector




Inserts*:



CEV398IU: 340 mm
CEV898IU: 240 mm

1-piece:




De Lagausie CEV998MU: 225 mm

Needle-holder

Straight tungsten jaws forceps:




CEV992B: 340 mm
Reinberg CEV992: 220 mm



CEV892: 220 mm

Curved tungsten jaws forcep:




CEV992C: 340 mm
CEV992D: 220 mm

*Inserts compatible with every rotative handles available on page 8 & 9

Scissors with serrated blades

Curved:



Inserts*:



CEV393IU: 330 mm
CEV893IU: 230 mm

Straight:



Inserts*:



CEV394IU: 330 mm
CEV894IU: 210 mm

1-piece:



CEV394BU: 350 mm
Montupet CEV894MU: 220 mm



De Lagausie CEV996MU: 220 mm

Integra® MicroFrance®

3.5 mm Laparoscopic Instruments

Rotative handles

Versatility: Interchangeable handles with all 3.5 mm and multiple 5 mm inserts*

Ergonomic handles



CEV720BRU: large grip - with declutching ratchet



CEV720RU: fine grip - with declutching ratchet



CEV7285RU: fine grip - with ratchet



CEV728BR5U: large grip - with declutching ratchet

Axial handles



CEV6385RU: no ratchet



CEV6285R: with ratchet (3 teeth)



CEV628P5: with ratchet (5 teeth)



CEV728T5: with external ratchet

* For more information, refer to laparoscopic catalogue #0335697

Angled handles



CEV1019T5BU: Large handle - with external ratchet - 2 mobile arms



CEV10195CU: with internal ratchet - 1 mobile arm



CEV1039T5TBU: large handle - no ratchet - 2 mobile arms



CEV10195RU: with external ratchet - 1 mobile arm



CEV10395DU: no ratchet - 1 mobile arm

Integra® MicroFrance®

3.5 mm Laparoscopic Instruments


Bipolar* Forceps – 1 piece instrument

High grasping capabilities and ergonomics, especially thanks to the « K » handle


Serrated jaws



Spatula serrated jaws



Fenestrated grasping jaws



* Compatible with the bipolar cables listed on page 15

Mini-Port Instruments

A large offer of Mini-Port Instruments completing the forceps product range

Monopolar hook

Round tip:



CEV390B: 330 mm
Montupet CEV890: 225 mm

90°:



CEV390C: 355 mm
CEV890C: 90° 225mm Ø 3,5mm

Suction irrigation tubes

Eyes suction irrigation tube, blunt, graduated 50/100/150mm:



Gayet CEV103G3: 295 mm
CEV8103G3: 200 mm

HANDLE



CEV103GH



Gayet CEV103BG

Eyes suction irrigation tube with blunt tip:




CEV896: 220 mm

Integra® MicroFrance®


3.5 mm Laparoscopic Instruments

Palpation probe




CEV203: Round probe 330mm

Knot pusher



Iannelli CEV342: 340 mm



Iannelli CEV342E: 340 mm


Cadière Forceps

Percutaneous fenestrated forceps with sharp jaws




Cadière CEV350: 225 mm

Retractable knife



CEV898: 190 mm



CEV898L: knife blade, (pack of 3)

Reusable trocars*


Reduced trocar diameter to limitate the scare size

Stainless steel family

Trocar tip:

Sharp:
50 mm: **CEV8143A3**
80 mm: **CEV143A3**


Blunt:
50 mm: **CEV8143A4**
80 mm: **CEV143A4**



Trocar head:

with silicone valve: **CEV127T3**

with magnetic valve: **CEV177A3**



Sleeve:

with stopcock:
50 mm: **CEV8143A7**
80 mm: **CEV143A7**

without stopcock:
50 mm: **CEV8143B7**
80 mm: **CEV143B7**

*For replacement accessories, please refer to the accessories listed in page 15


Integra® MicroFrance®

3.5 mm Laparoscopic Instruments


Plastic family

Trocar tip:

Sharp:
45 mm: **Montupet CEV8826P**
70 mm: **CEV326P**



Blunt:
45 mm: **CEV8826M**
70 mm: **CEV326M**



Trocar body:

smooth sleeve - no stopcock
45 mm: **Montupet CEV8826L**
70 mm: **CEV326L**



screw sleeve - no stopcock
45 mm: **Montupet CEV8826V**
70 mm: **CEV326V**



Accessories

Adapter cleaning luer-lock for tube insert



CEV39L
Adapter cleaning luer-lock for tube insert
(pack of 3)



Monopolar cables - 3.5m

Reference	Description
CP391E	For Erbe ACC ICC/VIO Storz Wolf Generators
CP391M	For Aesculap Berchtold Martin Erbe T Serie Lamidey Sutter BM-780 II Generators
CP391V	For Bowa Bovie Conmed Valleylab Generators

Bipolar cables – 4.5 m

Reference	Description
CP393-2	For Bowa Emc Lamidey Valleylab Generators
CP393-3	For Erbe Storz Wolf Generators
CP393-4	For Aesculap Berchtold Martin Sutter BM 780 II Generators

Trocar Accessories

Reference	Description	Quantity
CEV179A5	Replacement magnetic valves for trocar's head CEV177A3	6
CEV125P6	Replacement silicone valves for trocar's head	6
CEV103P6	Sealing Caps for 3.5mm Diameter Trocar's Head, Green	6

Sizes and scale are not contractual and are given for information purposes only. All the dimensions indicated in grey and in mm are for Working Length indications.

INDICATIONS FOR USE

All the medical devices showed in this brochure have indications for laparoscopy surgery.

CONTRAINDICATIONS

Do not use in indications other than the indications mentioned in the DFU or Labeling.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indications, contraindications, hazards, warnings, precautions, and instructions for use.

Products mentioned in this document are CE class I, IIa and IIb devices. Contact Integra should you need any additional information on devices classification. All the medical devices mentioned on this document are CE marked according to European council directive 93/42/EEC on medical devices and its relatives, unless specifically identified as "NOT CE MARKED".

For more information or to place an order, please contact:

Sales & Marketing EMEA

Integra LifeSciences Services (France) SAS
Immeuble Séquoia 2 ■ 97 allée Alexandre Borodine
Parc technologique de la Porte des Alpes
69800 Saint Priest ■ FRANCE
Phone: +33 (0)4 37 47 59 00 ■ Fax: +33 (0)4 37 47 59 30
integralife.eu

Manufacturers :



Integra MicroFrance SAS (France)
Le Pavillon ■ 03160 Saint Aubin le Monial ■ France



Integra Contact - Questions, Information & Product Ordering

International: +33 (0) 437 47 59 50 ■ +33 (0) 437 47 59 25 (Fax) ■ csemea@integralife.com
France: +33 (0) 437 47 59 10 ■ +33 (0) 437 47 59 29 (Fax) ■ custservfrance@integralife.com
United Kingdom: +44 (0) 1264 312 725 ■ +44 (0) 1264 312 821 (Fax) ■ custsvcs.uk@integralife.com
Italy: +39 (0)2 577 89 21 ■ +39 (0)2 575 113 71 (Fax) ■ custsvcitaly@integralife.com
Benelux: +32 (0)2 257 4130 ■ +32 (0)2 253 2466 (Fax) ■ custsvcbenelux@integralife.com