



Endonasal surgery product line

Integra® MicroFrance®

De Gabory endonasal surgery set

An ergonomic and flexible solution

Integra® MicroFrance®


De Gabory Instruments

Endonasal surgery set

- Incision-free access to the entire frontal sinus using micro-grasping forceps

- + Movement freedom thanks to the ergonomic stylet-like handle

- + Access to the anterior wall with various angled inserts




- Intuitive assembly and disassembly of inserts with audible click

- + Jaws




- + Interchangeable inserts

MCEN68CDG




MCEN68BDG



- + Stylet-like handle


MCEN728DG



45°




MCEN68ADG



- + Irrigation port accessory MCEN39L
(supplied in set of 3)




■ Suction controlled using the handle with Luer-lock connector



Catalog Number	Description (WL: Working length)
MCEN103LDG	Rotating connector handle for suction
MCEN68ADG	Grasping forceps 45° insert 3 x 6 mm with angled tip 120 mm WL for MCEN728DG handle
MCEN68BDG	Grasping forceps insert 3 x 6 mm with angled tip 120 mm WL for MCEN728DG handle
MCEN68CDG	Grasping forceps insert 3 x 6 mm with straight tip 120 mm WL for MCEN728DG handle
MCEN728DG	Stylet handle for grasping forceps insert

■ **K-handle makes the micro-grasping forceps easy to manipulate**

- + Provides a natural gripping action
- + Suited to all hand sizes
- + Makes it easier to pass instruments between operators



Grasping forceps	Catalog Number	WL 110 mm K-handle		
	MCEN33DG MCEN34DG	fenestrated	straight 45°	3 x 6 mm 3 x 6 mm
	MCEN62T1DG MCEN62A1DG MCEN62B1DG	fenestrated	straight 45° 90°	3 x 8 mm 3 x 8 mm 3 x 8 mm
	MCEN62T2DG MCEN62A2DG MCEN62B2DG	fenestrated	straight 45° 90°	5 x 11 mm 5 x 11 mm 5 x 11 mm
	MCEN90DG MCEN91DG MCEN91BDG	fenestrated	straight 45° 45°	3 x 7 mm 3 x 7 mm 3 x 10 mm
	MCEN4DG MCEN5DG MCEN5BDG		straight 45° 45°	3 x 6 mm 3 x 6 mm 3 x 12 mm

Indications: All the medical devices listed in this document are indicated for nasal surgery.

Contraindications: These instruments must only be used for the intended indications listed in the IFU or labeling.

Additional information for customers in Europe, the Middle East and Africa:

Caution: Federal Law restricts this device to sale by or on the order of a physician.

All the products listed in this document are CE-marked in accordance with the EU Medical Device Directive 93/42/EEC and its relatives, except for products specifically identified as "Not CE-marked".

All the text and diagrams in this catalog are illustrative and for information only. None of the information in this catalog is intended to act as a substitute for the instructions included in the packaging of every Integra product. Always refer to the appropriate IFU for detailed clinical instructions. For additional technical or clinical information on Integra products, please contact your Integra sales representative.

Carefully read the product's IFU. The dimensions listed in this document are not contractual and are provided for information only.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

Class I and IIa products. For more information on the product class, please contact Integra Customer Service.

Last revision date: 07/2017

For more information, please contact customer service:

Europe, Middle East, Africa

France / International : +33 (0) 470 67 98 00 ■ +33 (0) 470 67 98 20 (Fax) ■ microfrance@integralife.com
 United Kingdom : +44 (0) 1264 312 725 ■ +44 (0) 1264 312 821 (Fax) ■ custsvckuk@integralife.com
 Germany : +49 (0) 2102 5535 6200 ■ +49 (0) 2102 5536 636 (Fax) ■ custsvcgermany@integralife.com
 Benelux : +32 (0)2 257 4130 ■ +32 (0)2 253 2466 (Fax) ■ custsvcbenelux@integralife.com
 Switzerland : +41 (0)22 721 23 00 ■ +41 (0)22 721 23 99 (Fax) ■ custsvcsuisse@integralife.com
 Italy : +39 (0) 2 577 89 21 1 ■ +39 (0) 2 575 113 71 (Fax) ■ custsvcitaly@integralife.com
integralife.eu/contact

Integra and the Integra logo are registered trademarks of Integra LifeSciences in the United States and/or other countries

©2017 Integra LifeSciences Corporation. All rights reserved.

0393882-1-EN

DOCUMENT FOR USE IN EUROPE, MIDDLE EAST and AFRICA ONLY.

Manufacturer:



Integra MicroFrance SAS (France)
 Le Pavillon ■ 03160 Saint Aubin le Monial ■ France