

Integra® MicroFrance®

Sataloff Indirect Laryngeal
Instrument Set



Integra® MicroFrance®

Sataloff Indirect Laryngeal Instrument Set

Designed in conjunction with Robert T. Sataloff, MD

The Sataloff Indirect Laryngeal Instrument Set is used for office-based procedures. This set includes dilators, cupped forceps, scissors, heart-shaped forceps, alligator forceps, and suctions that are curved to help facilitate access to the laryngeal region with the patient in a sitting position. The smaller instrument tray has been designed for use with office-oriented sterilization equipment.


Sataloff Indirect Laryngeal Set

Reference	Description	Qty.
MCLS62**	Sataloff Indirect Laryngeal Instrument Set	1
Set Includes:		
MCLS44	Laryngeal Dilator, Size 0	1
MCLS45	Laryngeal Dilator, Size 1	1
MCLS46	Laryngeal Dilator, Size 2	1
MCLS47	Laryngeal Dilator, Size 3	1
MCLS48	Modular Handle	1
MCLS49	Indirect Cupped Forceps, Straight	1
MCLS50	Indirect Cupped Forceps, Left	1
MCLS51	Indirect Cupped Forceps, Right	1
MCLS52	Indirect Scissors, Straight	1
MCLS53	Indirect Scissors, Right	1
1650030	Oro-Tracheal Injector with Curved Cannula, 22cm	1
MCLS54	Indirect Scissors, Left	1
MCLS55	Indirect Heart-Shaped Forceps, Left	1
MCLS56	Indirect Heart-Shaped Forceps, Right	1
MCLS57	Indirect Alligator Forceps, Straight	1
MCLS58	Indirect Alligator Forceps, Right	1
MCLS59	Indirect Alligator Forceps, Left	1
MCLS60	Indirect Suction Tube, 5FR	1
MCLS61	Indirect Suction Tube, 7FR	1
MCLS63	Vascular Knife	3
MCLS64	Sharp Knife	3
MCLS65	Sickle Knife	3
MCLS66	Flap Knife	3
MCLS62SC	Instrument Tray, Stainless Steel	1

**Not for sale outside US



MCLS62



The MicroFrance® brand is a premier line of ENT instruments, designed with world-renowned ENT surgeons to facilitate new techniques and improve patient outcomes.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.


- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:


United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ■ 888-980-7742 fax
International +1 609-936-5400 ■ +1 609-750-4259 fax
integralife.com

Manufacturer:

 **Integra MicroFrance SAS**
Le Pavillon, 03160
Saint-Aubin-le-Monial ■ France

Manufacturer:

 **Integra LifeSciences Production Corporation**
11 Cabot Boulevard
Mansfield, MA 02048 USA