

Integra®

Table Mounted
Speed-Tract® Retractor System


Limit uncertainty when
speed is of the essence



INTEGRA®
LIMIT UNCERTAINTY

Integra® Speed-Tract®

Simple two piece design.




Speed-Tract® Sterile Field Post



Speed-Tract® Support Arm

Simple Setup for Large Patients




Quick Locking Mechanism




The Integra Speed-Tract® has a locking mechanism which locks and unlocks as simply as a lightswitch.

Example of Competitive Product


Three piece design.



Patient Setup



Threaded Connectors



Integra Speed-Tract® Ratchets

Infinite retractor positioning capability in every plane



#11646 Tilt Ratchet Clamp

- Standard design with tilt



#11560 Integrated Tilt Ratchet with Jaw

- Incorporates a lower jaw to keep the clamp on the ring




#11490 Ratchetball™

- Ball and socket joint allows infinite retractor positioning

Integra Speed-Tract® Split Frame

The Split Frame is a one-piece adjustable frame that allows for low profile or elevated retraction.



Integra Speed-Tract® Component Parts



Integra Speed-Tract® Component Parts List

Interchangeable with all marketed Ring Retractor Systems.

Framing Components

Catalog Number	Description
11184	Support Arm
11321	Sterile Field Post, w/Coupler, 18"
11352	Sterile Field Post, w/Coupler, 24"
11275	Coupler, 1" x 1"
11826	1" DIA Vertical Extension Bar

Rings

Catalog Number	Description
11631	Extra Small Oval Ring, 6" x 7 ³ / ₄ "
11632	Small Oval Ring, 8 ¹ / ₂ " x 12"
11481	Medium Oval Ring, 10 ¹ / ₂ " x 15"
11482	Large Oval Ring, 12 ¹ / ₂ " x 18 ¹ / ₂ "
11633	Small Round Ring, 8 ¹ / ₂ "
11480	Medium Round Ring, 10 ¹ / ₂ "
11634	Large Round Ring, 12 ¹ / ₂ "
11635	Small Segmented Ring, 8 ¹ / ₂ " x 14"
11875	Medium Segmented Ring, 10 ¹ / ₂ " x 19 ¹ / ₂ "
11637	Round Segmented Ring, 8 ¹ / ₂ "
11966	Segmented Perineal Ring
11968	Rigid Perineal Ring
11848	Split Frame, Oval
11849	Split Frame, Round

Clamps:

Catalog Number	Description
11910	Standard Blade Clamp
11646	Tilting Ratchet Clamp, w/o lip
11560	Tilting Ratchet Clamp, with lip
11490	Ratchet Ball Clamp

Balfour Retractors:

Catalog Number	Description
11902	Balfour Fixed Retractor, SQ Shaft, 2 ¹ / ₄ " x 1 ¹ / ₂ "
11510	Balfour Fixed Retractor, SQ Shaft, 2 ³ / ₄ " x 3"
11942	Balfour Fixed Retractor, SQ Shaft, 3 ¹ / ₂ " x 2 ³ / ₄ "
11440	Balfour Swivel Retractor, SQ Shaft, 2 ³ / ₄ " x 3"

Deaver Retractor:

Catalog Number	Description
11904	Deaver Fixed Retractor, SQ Shaft, 3 ⁴ / ₈ " x 3"

Gelpi Retractor:

Catalog Number	Description
11908	Gelpi Retractor, SQ Shaft

Harrington Retractors:

Catalog Number	Description
11511	Harrington Fixed Retractor, SQ Shaft, 2 ¹ / ₂ " x 7"
11952	Harrington Fixed Retractor, SQ Shaft, 3 ¹ / ₂ " x 10"
11441	Harrington Swivel Retractor, SQ Shaft, 2 ¹ / ₂ " x 7"

Kelly Retractors:

Catalog Number	Description
11645	Kelly Fixed Retractor, SQ Shaft, 1" x 2 ¹ / ₂ "
11644	Kelly Fixed Retractor, SQ Shaft, 1" x 1 ¹ / ₂ "
11640	Kelly Fixed Retractor, SQ Shaft, 2" x 2 ¹ / ₂ "
11512	Kelly Fixed Retractor, SQ Shaft, 2" x 3"
11513	Kelly Fixed Retractor, SQ Shaft, 2" x 4"
11514	Kelly Fixed Retractor, SQ Shaft, 2" x 5"
11515	Kelly Fixed Retractor, SQ Shaft, 2" x 6"
11641	Kelly Fixed Retractor, SQ Shaft, 2" x 8"
11442	Kelly Swivel Retractor, SQ Shaft, 2" x 3"
11443	Kelly Swivel Retractor, SQ Shaft, 2" x 4"
11444	Kelly Swivel Retractor, SQ Shaft, 2" x 5"
11445	Kelly Swivel Retractor, SQ Shaft, 2" x 6"

Malleable Retractors:

Catalog Number	Description
11938	Malleable Fixed Retractor, SQ Shaft, 1/2" x 4"
11940	Malleable Fixed Retractor, SQ Shaft, 1" x 5"
11516	Malleable Fixed Retractor, SQ Shaft, 1 1/2" x 6"
11944	Malleable Fixed Retractor, SQ Shaft, 2" x 10"
11517	Malleable Fixed Retractor, SQ Shaft, 2" x 6 1/2"
11950	Malleable Fixed Retractor, SQ Shaft, 2" x 8"
11946	Malleable Fixed Retractor, SQ Shaft, 3" x 10"
11518	Malleable Fixed Retractor, SQ Shaft, 3" x 6 1/2"
11948	Malleable Fixed Retractor, SQ Shaft, 4" x 8"
11446	Malleable Swivel Retractor, SQ Shaft, 1 1/2" x 6"
11447	Malleable Swivel Retractor, SQ Shaft, 2" x 6 1/2"
11448	Malleable Swivel Retractor, SQ Shaft, 3" x 6 1/2"
11820	Notched Malleable Fixed Retractor, SQ Shaft, 2" x 6"
11822	Notched Malleable Fixed Retractor, SQ Shaft, 2" x 8"
11821	Notched Malleable Fixed Retractor, SQ Shaft, 2 1/2" x 7"

Renal Retractor:

Catalog Number	Description
11956	Renal Fixed Retractor, SQ Shaft, 1" x 6"

Richardson Retractors:

Catalog Number	Description
11954	Richardson Fixed Retractor, SQ Shaft, 3/4" x 1"
11936	Richardson Fixed Retractor, SQ Shaft, 3/4" x 2"

Vaginal Retractors:

Catalog Number	Description
11958	Vaginal Lateral Fixed Retractor, SQ Shaft, 1" x 3"
11962	Vaginal Deaver Fixed Retractor, Large, SQ Shaft, 1" x 5"
11960	Vaginal Deaver Fixed Retractor, Small, SQ Shaft, 1" x 3 1/2"
11964	Vaginal Posterior Fixed Retractor, SQ Shaft, 1" x 4 1/2"

Aluminum Blades:

Catalog Number	Description
11300R	Aluminum, Deaver ABH Blade, 2 ³ / ₄ " x 4 ¹ / ₂ "
11301R	Aluminum, Malleable ABH Blade, 2" x 6 ¹ / ₂ "
11302R	Aluminum, Mayo ABH Blade, 2 ³ / ₄ " x 2"
11303R	Aluminum, Mayo ABH Blade, 3 ¹ / ₂ " x 2"
11304R	Aluminum, Renal ABH Blade, 1" x 4"
11305R	Aluminum, Renal ABH Blade, 1" x 6"
11306	Aluminum, Renal ABH Blade, 1" x 8"
11308	Aluminum, Splanchnic ABH Blade, 2" x 4 ¹ / ₂ "
11309	Aluminum, Splanchnic ABH Blade, 2" x 6"
11863	Aluminum, Renal ABH Blade, 1" x 3"
11864	Aluminum, Renal ABH Blade, 1" x 4"
11865	Aluminum, Renal ABH Blade, 1" x 5"
11866	Aluminum, Renal ABH Blade, 1" x 6"
11867	Aluminum, Renal ABH Blade, 1" x 7"
11868	Aluminum, Renal ABH Blade, 1" x 8"

Blade Holder for ABH Style Blades:

Catalog Number	Description
11340	Square Blade Holder

*ANSI/AAMI ST77:2006 4.3.5 "Maximum Weight"

AORN "Recommended practices for the selection and use of packaging systems for sterilization" Recommended Practice I, Section 4

Integra®

Speed-Tract Retractor System

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-367-8657 ■ 651-287-4400 fax

Outside USA 651-287-4300 ■ 651-287-4400 fax

integralife.com

