

CUSA® Clarity System


Designed for Nurses





Simplified set up¹
Everything is ready when the surgeon needs it

Intuitive interface
Quick responsiveness in the OR


QUICK & SIMPLE
CUSA® Clarity



- + Streamlined console setup thanks to guided step-by-step instructions available right on the touchscreen
- + Faster irrigation priming¹
- + One-step tubing assembly with CUSA Quick Connect™ cartridge and tubing set



- + Hexacog handpiece design prevents overtorquing to ensure reliability and optimal performance



- + Optional proportional footswitch mode enables surgeon control of amplitude
- + Enhanced troubleshooting audible alert sounds and text displays on screen
- + On-screen help provides real-time assistance on demand



- + Simple fine-tuning of the settings
Amplitude
TissueSelect®
Aspiration
Irrigation

¹ Compared to CUSA® Excel+ system

Ordering information

Capital Equipment

Catalog No.	Description	Sterile	Quantity
C7000	CUSA® Clarity Console with Footswitch	No	Each 1
C7001	Cart	No	Each 1
C7036	36 kHz Handpiece with Cleaning Brush	No	Each 1

Tips

Catalog No.	Description	Tip length	Tip ID*	Sterile	Quantity
Conventional Single Use Tips					
C7411S	36 kHz Curved Extended MicroTip™	12,15 cm (4.78 in)	1,6 mm (0.06 in)	Yes	Box 5
C7411SEA					Each 1
C7412S	36 kHz Standard Tip Short	4,57 cm (1.80 in)	2,0 mm (0.08 in)	Yes	Box 5
C7412SEA					Each 1
C7414S	36 kHz Curved Extended Standard Tip	11,44 cm (4.50 in)	2,0 mm (0.08 in)	Yes	Box 5
C7414SEA					Each 1
Specialty Single Use Tips					
C7415S	36 kHz MicroTip™ Plus	19,27 cm (7.58 in)	1,6 mm (0.06 in)	Yes	Box 5
C7415SEA					Each 1
C7417S	36 kHz ShearTip™	12,05 cm (4.74 in)	1,6 mm (0.06 in)	Yes	Box 5
C7417SEA					Each 1

Consumable

Catalog No.	Description	Sterile	Quantity
C7300	CUSA Quick Connect™ Tubing	Yes	Box 5
C7300EA		Yes	Each 1
C7602	36 kHz Torque Wrench	Yes	Box 5
C7602EA		Yes	Each 1

Accessories

Catalog No.	Description	Sterile	Quantity
C7636	36 kHz Reusable Torque Base	No	Each 1
C7136	36 kHz Sterilization Tray	No	Each 1
C7005	Contamination Guards	No	Box 2
C7500	Cleaning Brush	No	Each 1
C7550	Specimen Trap	No	Box 12
C7002	Footswitch	No	Each 1

Indications: The CUSA Clarity Ultrasonic Surgical Aspirator system is indicated for use in these surgical procedures where fragmentation, emulsification and aspiration of soft tissue is desirable: Neurosurgery, Gastrointestinal and affiliated organ surgery, Urological surgery, Plastic and reconstructive surgery, General surgery, Orthopedic surgery, Gynecological surgery, Thoracic surgery, Laparoscopic surgery and Thoracoscopic surgery.

Contraindications: This ultrasonic surgical aspirator device is not indicated for and should not be used for the removal of uterine fibroids.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Consult product labels and inserts for any indications, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

Integra Canada ULC

2590 Bristol Circle, Unit 1

Oakville, Ontario ■ Canada

L6H6Z7

1-800-263-6693 (phone) ■ 905-632-7938 (fax)

custserv.canada@integralife.com



Manufacturer:



Integra LifeSciences (Ireland) Ltd
IDA Business and Technology Park
Sragh, Tullamore, County Offaly ■ IRELAND