


Catheters



Peritoneal Catheters, Accu-Flo Barium Distal Slit Valve


The slit valve in the Codman® Neuro distal valve catheters are 90° apart and are precisely cut to achieve more consistent closing pressures. The dots adjacent to the slits indicate the pressure.

Catalog#	Closing Pressure (Measured in cm of water)	Length	Description
82-1301	2-5 cm H ₂ O Low	91 cm, 1.27 mm ID	• 4 pressure-regulating slits
82-1302	5-9 cm H ₂ O Medium	x 2.54 mm OD	• Tantalum marks and tip
82-1303	9-14 cm H ₂ O High		


Ventricular Catheters

All ventricular catheters are packaged sterile and are made from silicone tubing.

Catalog#	Description	Includes:
82-1201	Accu-Flo® Barium Ventricular Catheter	<ul style="list-style-type: none"> • 15 cm, 1.3 mm ID x 2.5 mm OD • X-ray detectable throughout. Rounded X-ray detectable tantalum tip • X-ray detectable tantalum markings, 5 cm and 10 cm. With 40 inlet holes, four rows of ten holes, 90° apart • Stainless steel stylet
82-1203	Accu-Flo Clear Ventricular Catheter	<ul style="list-style-type: none"> • 30 cm, 1.3 mm ID x 2.5 mm OD • Rounded X-ray detectable tantalum tip • X-ray detectable tantalum mark at 15 cm • With 40 inlet holes, four rows of ten holes, 90° apart
82-1221	Accu-Flo Dotted Ventricular Catheter	<ul style="list-style-type: none"> • 15 cm, 1.3 mm ID x 2.5 mm OD • Rounded X-ray detectable tantalum tip • Clear silicone tubing with X-ray detectable tantalum markings every 1 cm • With 40 inlet holes, four rows of ten holes, 90° apart • Stainless steel stylet




Accu-Flo Barium Right Angle Ventricular Catheter



Catalog#	Length	Description	Includes:
82-1231	5 cm	Accu-Flo Right Angle Ventricular Catheter	• 4 rows of 10 holes
82-1232	7 cm		• 1.27 mm ID x 2.54 mm OD
82-1233	9 cm		• Tantalum tip
82-1234	11 cm		• Stainless steel stylet

Holter® Barium Ventricular Catheter



Catalog#	Length	Number of Perforations	Description	Includes:
82-1652	3 cm	3	Ventricular Catheter, with a Right Angle	• 1.5 mm ID x 3.1 mm OD
82-1654	4 cm	3		• All ventricular tips have multiple perforations
82-1658	6 cm	4		• Stainless steel stylet
82-1660	7 cm	5		
82-1662	8 cm	5		

Catalog#	Description	Includes:
82-1650	Ventricular Catheter, straight	• 15 cm, 1.5 mm ID x 3.1 mm OD • Stainless steel stylet






Codman Ventricular Catheter

Catalog#	Description
82-3041	Codman MEDOS® Ventricular Catheter <ul style="list-style-type: none"> • Silicone, impregnated with barium • Right-angled adaptor, plastic • 14 cm, 1.4 mm ID x 2.7 mm OD • 3 rows of 8 holes (.5 mm hole diameter) • Stainless steel stylet

Atrial Catheters

Holter Distal Atrial Catheters, sterile, barium impregnated with two lumen configurations to meet various shunting needs. Large lumen styles have an ID of 1.2 mm and an OD of 2.5 mm. The small, or stepped down portion lumen have an ID of .8 mm and an OD of 1.3 mm.







Catalog#	Description
82-1670	Holter Distal Atrial Catheter, Type A, barium,sterile. Large lumen catheter. <ul style="list-style-type: none"> • Total length 45 cm • Two sterile 82-1694 Holter Type A Connectors, stainless steel



Catalog#	Description
82-1676	Holter Distal Atrial Catheter, Type E, barium, sterile <ul style="list-style-type: none"> • Small catheter length 24 cm • Large catheter length 24 cm • Total length 48 cm • One sterile 82-1694 Holter Type A Connector, stainless steel



Catalog#	Description
82-3044	Codman Atrial Catheter <ul style="list-style-type: none"> • 46 cm, 1 mm ID x 2.2 mm OD



Holter Distal Peritoneal Catheters

Holter Distal Peritoneal Catheters, Barium, Sterile. Barium impregnated with 1.2 mm ID x 2.5 mm OD. Supplied with two sterile Holter Type A Connectors (ref. 82-1694)

Catalog#	Description
82-1682	Holter Distal Peritoneal Catheter
	<ul style="list-style-type: none"> • Open end • Two sterile 82-1694 Holter Type A Connectors, stainless steel

Catalog#	Description
82-1684	Holter Distal Peritoneal Catheter
	<ul style="list-style-type: none"> • Closed tip with 4 side slits 90° apart • Tip is solid silicone and is more heavily opacified for positive X-ray identification • Two sterile 82-1694 Holter Type A Connectors, stainless steel

MEDOS Peritoneal Catheter

Catalog#	Description
82-3045	Codman Peritoneal Catheter
	<ul style="list-style-type: none"> • 120 cm, 1 mm ID x 2.2 mm OD • Barium impregnated silicone

Accu-Flo Peritoneal Catheters

Catalog#	Description
82-1380	Accu-Flo Barium Distal Open-end Catheter
	<ul style="list-style-type: none"> • Length 91 cm • Made of barium impregnated silicone tubing • 1.3 mm ID x 2.5 mm OD X-ray detectable • Each catheter has 24 drainage holes, four rows of six holes, 90° apart • X-ray detectable tantalum markings are at the distal 10 cm, 20 cm, and 30 cm

Page intentionally left blank

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order, please contact:

United States, Canada, Asia, Pacific, Latin America

USA 800-654-2873 ■ 888-980-7742 fax

International +1 609-936-5400 ■ +1 609-750-4259 fax

integralife.com

