


# Codman®

Surgical Patties and Strips

---

**COMPARE THE QUALITY.  
EXPERIENCE THE DIFFERENCE.**



**Codman®**  
SPECIALTY SURGICAL

## The Benchmark in Quality


**Description** - Codman Surgical Patties and Strips provide the neurosurgeon with softness, gentleness, protection, absorbency, and convenience to meet the demands of delicate intracranial and spinal surgery.

**Absorbent** - Codman Patties and Strips absorb more than 7 times their own weight in less than a second.

**X-Ray Detectable** - Codman Patties and Strips have built-in x-ray-detectable markers.

**Protection** - Codman Patties and Strips are made of a proprietary Cottonoid® material of compressed rayon woven by a special process.

**Easy to Remove** - Sutureweld strong bonds hold the indicator string to the patty while the surface that lies next to tissue is completely smooth and free of stitches.



**Indications** - The Surgical Patties and Surgical Strips are indicated for use in protection of tissue, including brain and other tissues of the central nervous system, during surgery.

Please see Instructions for Use for additional warnings.

## Ordering Information

### Codman® Surgical Patties

(Sterile, X-ray Detectable)

10 patties/package, 20 packages/dispenser box

Reference	Size	
80-1396	5/16" Round	(7 mm) Round
80-1397	1/4" x 1 1/2"	(7 mm x 38 mm)
80-1398	1/4" x 3"	(7 mm x 76 mm)
80-1399	1/4" x 1/4"	(7 mm x 7 mm)
80-1400	1/2" x 1/2"	(13 mm x 13 mm)
80-1401	3/4" x 3/4"	(19 mm x 19 mm)
80-1402	1/2" x 1"	(13 mm x 25 mm)
80-1403	1" x 1"	(25 mm x 25 mm)
80-1404	1/2" x 1 1/2"	(13 mm x 38 mm)
80-1406	1/2" x 2"	(13 mm x 52 mm)
80-1407	1/2" x 3"	(13 mm x 76 mm)
80-1408	1" x 3"	(25 mm x 76 mm)
80-1409	3" x 3"	(76 mm x 76 mm)

### Codman® Surgical Strips

(Sterile, X-ray Detectable)

10 strips/package, 20 packages/dispenser box

Reference	Size	
80-1449	1/8" x 6"	(3 mm x 152 mm)
80-1450	1/4" x 6"	(6 mm x 152 mm)
80-1451	1/2" x 6"	(13 mm x 152 mm)
80-1452	3/4" x 6"	(19 mm x 152 mm)
80-1453	1" x 6"	(25 mm x 152 mm)
80-1454	1 1/2" x 6"	(38 mm x 152 mm)
80-1455	2" x 6"	(51 mm x 152 mm)
80-1456	3" x 6"	(76 mm x 152 mm)
80-1457	3 1/2" x 6"	(89 mm x 152 mm)

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region.

- Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality.
- Warning: Applicable laws restrict these products to sale by or on the order of a physician.
- Consult product labels and inserts for any indication, contraindications, hazards, warnings, precautions, and instructions for use.

For more information or to place an order please contact:

#### United States

USA 800-654-2873 ▪ 888-980-7742 fax

International +1 609-936-5400 ▪ +1 609-750-4259 fax

[integralife.com](http://integralife.com)

